

What is PRESEPT™ ?

PRESEPT™ is recommended for **disinfection of :**

- ✓ Floors, walls, work surfaces, cupboards
- ✓ Baby bottles and other utensils
- ✓ Soiled or infected linen
- ✓ Clean up of body fluids and blood spillage (PRESEPT™ Granules)
- ✓ Non critical items requiring disinfection (e.g. bedpans)



PRESEPT™ – Features & Benefits

Product features	Benefits
Contains NaDCC	Broad-spectrum activity More effective than hypochlorites
Stable solution at acidic pH	Better conditions for chlorine biocidal power
Unique “50 : 50” equilibrium properties	Extended biocidal activity – 24h reuse Better resistance to organic matter
Tablets and Granules forms	Convenient to use – Compact to store
Clear dilution instructions for Tablets	Multi-using purpose - Economical
Available in different strengths and packaging	Suitable for every hospital need

Why is PRESEPT™ more effective than hypochlorites?

Unique « 50 : 50 » equilibrium properties of PRESEPT™

- ✓ NaDCC dissolves in water to release HOCl - measured as available chlorine. Only 50% of the chlorine is actually "freed" to kill micro-organisms. The other 50% stay "combined". As "free available" chlorine is used up (i.e. killing microorganisms), "combined" chlorine is continuously released to restore the permanent 50 : 50 ratio.



Results : prolonged bactericidal effect and increased resistance to organic matter.

- ✗ Hypochlorites don't have such a feature, which makes the diluted solution less effective over time, even at the same buffered pH.

PRESEPT™ – Materials Compatibility

PRESEPT™ Tablets and Granules can be used to disinfect **hard surfaces or items** made from:

- Glassware
- Rubber
- Plastic
- Porcelain
- High-grade stainless steel
- Linen (not made from animals)