




3M™ Attest™ Mini Auto-reader 490M System* Quick Reference Guide

VH2O2 BI Preparation and Use



1 Inspect
Remove BI from pouch and reseal. Verify media ampoule is intact and indicator stripes on cap are blue.



2 Pouch
Place BI and chemical indicator (CI) in peel pouch indicated for use in vaporized hydrogen peroxide (VH2O2) sterilization processes.



3 Process
Seal pouch and place in sterilizer with white side facing up. Process load as recommended by sterilizer manufacturer.



4 Identify
After completion of cycle, inspect and verify media ampoule is intact and indicator stripes on cap have turned towards pink. Identify the BI by writing on the label.

STEAM and VH2O2 BI Activation and Incubation



1 Crush
Squeeze BI in activator to close cap and crush media ampoule.



2 Flick
Remove BI and flick wrist to ensure media flows to growth chamber.



3 Inspect
Verify media has transferred to BI growth chamber at the vial bottom.




4 Incubate
Place activated BI in any well. Display shows remaining minutes of incubation or test results.
DO NOT REMOVE BI until the (+) or (-) symbol indicates the test is complete.

3M™ Attest™ Mini Auto-reader 490M

Positive Control BI
For both Steam and VH2O2, each day a processed BI is incubated, activate and incubate an unprocessed BI as a positive control. The control BI should be from the same lot code as the processed BI.

Results After Incubation
Control BI: Positive (+) **Acceptable**
Processed BI: Negative (-) **Acceptable Sterilization**
Processed BI: Positive (+) **Failed Sterilization**

Remaining Incubation Minutes



- NOTES:**
- Activate and incubate 1295 BI within 1 hour of completion of the VH2O2 sterilization cycle.
 - Gloves are recommended when activating 1295 BI.
 - Allow 1491 and 1492V BIs to cool 10 minutes prior to activation.

Interpretation of 3M™ Sterilization Assurance Core Products

| STEAM Indicator | Unexposed | | Acceptable Exposed | |
|---------------------------------|---|--|--------------------|--|
| 3M™ Attest™ BI 1491 and 1492V | | | | |
| | After processing, indicator on cap turns light brown or darker. | | | |
| 3M™ Comply™ Indicator Tape 1322 | | | | |
| | After processing, indicator stripes turn dark-brown/black. | | | |
| 3M™ Comply™ Indicator Tape 1355 | | | | |
| | After processing, indicator stripes turn dark-brown/black. | | | |

| VH2O2 Indicator | Unexposed | Acceptable Exposed Color Range | | |
|-------------------------------------|---|--------------------------------|--|--|
| 3M™ Attest™ BI 1295 | | | | |
| | After processing, indicator stripes turn toward pink. | | | |
| 3M™ Comply™ Indicator Tape 1228 | | | | |
| | After processing, indicator stripes turn toward pink. | | | |
| 3M™ Comply™ Chemical Indicator 1248 | | | | |
| | After processing, the blue bar changes to the color of the color match or pink. | | | |

| STEAM Indicator | Unprocessed | Acceptable Processed PASS | | Unacceptable Processed FAIL | |
|---|-------------|--|--|--|--|
| 3M™ Attest™ Steam Chemical Integrator 1243 | | | | | |
| | | After processing, indicator bar crosses into the ACCEPT region. | | If the color bar is in the REJECT region or on the line, it is considered a fail. | |
| 3M™ Comply™ Bowie-Dick Plus Test Pack 00135LF | | | | | |
| 3M™ Comply™ Early Warning Test Sheet (above) | | | | | |
| 3M™ Comply™ Bowie-Dick Test Sheet (below) | | | | | |
| | | When indicator ink circle on Early Warning Test sheet is lighter than the color standard, this may indicate that sterilizer maintenance should be scheduled. | | | |
| | | After processing, a uniform color change to dark brown/black is a pass. | | Fail result is indicated by a non-uniform color development, with a lighter colored area, usually near the center. | |

