- Do not immerse animal fibres such as wool or silk due to strong proteolytic action.
 Avoid prolonged contact with stainless steel.
 Warning : The rate of chlorine generation is accelerated in acid conditions.

Disinfection of	Required concentration of available chlorine	DILUTION RATES			Additional instructions
		0.5g tablets	2.5g tablets	5.0g tablets	
Blood Spillage†	10,000 ppm	18 tablets in 0.5l water	7 tablets in 1l water	9 tablets in 2.5l water	Pour over blood Using gloves, wipe up with disinfectant saturated disposable cloth
Pipette jars †	2,500 ppm	9 tablets in 1l water	9 tablets in 5l water	9 tablets in 10l water	Drop tablets into water filled pipette jar. Discard daily
General laboratory/ environmental use†	1,000 ppm	4 tablets in 1l water	4 tablets in 5l water	3.5 tablets in 10l water	Wipe down surfaces with disposable disinfectant saturated cloth
Baby bottles/teats, stainless steel utensils, porcelain, glassware, rubber and plastic tubing	140 ppm	1 tablet in 2l water	1 tablet in 10l water	1 tablet in 20l water	Immerse for 1 hr
Eating utensils and crockery	140 ppm	1 tablet in 2l water	1 tablet in 10l water	1 tablet in 20l water	Rinse for 1-2 mins
Soiled & Infected Linen	140 ppm	1 tablet in 2l water	1 tablet in 10l water	1 tablet in 20l water	Immerse for 1 hr prior to washing
Work surfaces, cupboards, floors, etc	140 ppm	1 tablet in 2l water	1 tablet in 10l water	1 tablet in 20l water	Wash down
Dishcloths, mops, etc	60 ppm	1 tablet in 4.6l (1 gal) water	1 tablet in 23l (5 gals) water	1 tablet in 46l (10 gals) water	Soak to bleach, clean and deodorise

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