

# Multifunction Wheelchair Scale with two ramps

H351-00-4 (Class III); H350-00-4 (Non-Class III)

Equipped with an integrated seat, **Shekel's** Electronic Multifunction Wheelchair Scale extends the concept of multifunctionality.

The accessible scale allows patients to be weighed while sitting on the safe and stable seat, using their own wheelchair, or even while standing supported by the robust handrail.

## Features

- Movement compensation technology
- BMI function
- Tare/Zero/Hold/Print/kg/lb
- Scroll Down/Up
- 5 digit ¾" LCD/Positive FSTN
- Automatic Reweigh
- Updatable firmware
- Ports: RS232 and USB
- Optional: Bluetooth

## Power

- 230v/120v input, DC 9V output
- Auto power off
- 6 X 1.5AA batteries
- Rechargeable battery (Class III only)
- Low battery indicator

## Technical Data

Class III	
Capacity	300kg / 660lb
Graduation	0.1kg / 0.2lb
Total Weight	43.2kg / 95.2lb
Product Size	1250 x 740 x 1250 mm 49.2" x 37" x 47.2"
Platform Size	800 x 780 x 70 mm 31.5" x 30.7" x 2.7"
Non-Class III	
Capacity	360kg / 800lb
Graduation	0.1kg / 0.2lb

Class III EMR READY



- In spite of its robustness and versatility, the scale can be easily folded



**SHEKEL** All the weigh!

**healthy medical**  
www.hmed.com.hk  
Sales Hotline: 8336 3333