



# Chair Scale

H551-00-1 (Class III); H550-00-1 (Non-Class III)

The Chair Scale with its patented design maintains **Shekel's** well-known reputation for reliability and uncompromising precision. With a few smart changes, **Shekel** has created the ideal ergonomic patient product. All of a patient's care needs, including safety, comfort, caregiver's use and maintenance, are met.

The compact package design ensures the most cost efficient shipping method.

## Features

- Movement compensation technology
- BMI function
- Tare/Zero/Hold/Print/kg/lb
- Scroll Down/Up
- 5 digit ¾" LCD/Positive FSTN
- Automatic Reweigh
- Updatable firmware
- Ports: RS232 and USB
- Optional: Bluetooth

## Power

- 230v/120v input, DC 9V output
- Auto power off
- 6 X 1.5AA batteries
- Rechargeable battery (Class III only)
- Low battery indicator

## Technical Data

Class III & Non-Class III	
Capacity	300kg / 660lb
Graduation	0.1kg / 0.2lb
Total Weight	18kg / 40lb
Product Size	630 x 850 x 965 mm 24.8" x 33.5" x 38"
Seat Size	430 x 420 x 80 mm 17" x 16.5" x 3"

Class III **EMR**  
READY

## Ergonomic patient-oriented design:

The armrests can be lifted separately, allowing the patient to slide from his chair to the chair scale while leaning on the prone armrest. The footrest can be folded as well for easier access.

