3M[™] Littmann[®] Stethoscopes

Introducing the next generation of a clinical icon.



Open bell versatility.

The pediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim.

Introducing the 3M[™] Littmann[®] Classic III[™] Stethoscope

Right out of the box, it lives up to its name.

The Littmann Classic III stethoscope is the latest version of the stethoscope that helps millions of medical professionals achieve their best. It brings new design, materials, and technology to the Littmann stethoscope series trusted by clinicians around the world. The new Classic III stethoscope offers you:



High acoustic sensitivity for performing general physical assessments.



Tunable diaphragms on both the adult and pediatric sides of the chestpiece.



Single-piece diaphragm that's easy to attach and easy to clean because the surface is smooth without crevices.



Next-generation tubing with improved resistance to skin oils and alcohol for longer tubing life. Less likely to pick up stains.

Using a Classic III stethoscope demonstrates what's important to you as a medical professional: quality, achievement, and providing the best possible care. Every patient. Every day.





Take a closer look at the 3M[™] Littmann[®] Classic III[™] Stethoscope

Features

Tunable technology is a 3M invention that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece.

Soft-sealing eartips give you a comfortable fit and an excellent acoustic seal; they snap tightly onto the eartubes for safety.

Open bell stays clear of dirt and debris by keeping it covered with the small diaphragm.

Adjustable headset lets you easily set the tension for head size and comfort just by squeezing or pulling apart the eartubes.

No natural rubber latex or phthalate plasticizers are used in the tubing or any other component—to help protect human health and the environment.

Made in the USA for consistent quality. We apply rigorous process and quality controls to every Littmann stethoscope. You'll always get a durable, dependable, clinical tool that will serve you for years to come.

Classic III stethoscope specifications 69cm / 27.2 inches Length Weight 150g Acoustic 7 performance Warranty 5 years CHESTPIECE Construction Two-sided, adult and pediatric Material Stainless steel Diaphragm type Tunable, single-piece Adult: 4.3cm / 1.7 in **Diaphragm diameter** Pediatric: 3.3cm / 1.3 in TUBING Single lumen Construction (single acoustic channel within tube) Not made with natural rubber latex

Wide diameter aerospace alloy /
anodized aluminum

or phthalate plasticizers

Colors and Finishes

We periodically introduce new tubing colors and chestpiece finish combinations for the Classic III stethoscope. Visit **littmann.com** for the latest product information.

Due to limitations of the printing process, the colors shown are approximations of the actual tubing colors.



Available tubing colors and finish combinations

TUBING COLOR	STANDARD FINISH	SPECIAL FINISH		
Black	5620	Black 5803, Smoke 5811		
Gray	5621			
Navy Blue	5622			
Caribbean Blue	5623	Rainbow 5807		
Hunter Green	5624			
Raspberry	5626	Rainbow 5806		
Burgundy	5627			
Orange	5629			
Ceil Blue	5630			
Rose Pink	5631*, 5639			
Pearl Pink	5633			
Chocolate		Copper 5809		
Dark Olive Green		Smoke 5812		
Plum	5831			
Lavender	5832			
Turquoise	5835			
Lemon-Lime	5839			
Lime Green	5829			
Emerald	5840			

*Available in the U.S. only

Material

HEADSET

Material



3M Health Care 3M Infection Prevention Division 2510 Conway Ave., St. Paul, MN 55144 USA

1-800-228-3957 | www.littmann.com

3M, Littmann, the Littmann Quality seal, the L, and Classic III are trademarks of 3M, used under license in Canada.

Please recycle. Printed in USA. © 3M 2015. All rights reserved. Made in USA. 70-2011-5600-0