

3M Science.
Applied to Life.™



Know now.

**Protect your patients
and your practice with
biological indicator
results in 24 minutes.**



3M™ Attest™ Mini Auto-reader 490M
for sterilization monitoring

Now available! The 3M™ Attest™ Mini Auto-reader.

In-office monitoring (spore testing) can help reduce the risk of patient exposure to contaminated dental instruments. Sized and priced for office use, the 3M™ Attest™ Mini Auto-reader 490M lets you bring rapid sterilization monitoring to your practice — because your patients deserve the highest standard of care.

- In-house sterilization monitoring results in just 24 minutes
- Know results before you release instruments for use — offering your patients a high standard of care
- Designed to be sized and priced for your practice
- Same expert 3M sterilization assurance technology as used in hospitals

Your partner in patient safety

- Easy to read display
- Small footprint for space-constrained areas
- Connect to a PC for digital documentation of biological indicator results



Catalog No.	Product	Other Notes	Packaging
490M	3M™ Attest™ Mini Auto-reader	Dimensions (W x H x D): 15,5cm x 8,25cm x 5,25cm (6.1 in x 3.2 in x 2.1 in)	1 unit/box
1491*	3M™ Attest™ Super Rapid Readout Biological Indicator (Blue cap) for Steam	Gravity: 132°C (270°F) 3 min. and 10 min. 135°C (275°F) 3 min. and 10 min.	50 each/box (4 boxes/case)
1492V*	3M™ Attest™ Super Rapid Readout Biological Indicator (Brown cap) for Steam	Dynamic-air-removal (Pre-vacuum and SFPP): 132°C (270°F) 3 min. and 4 min. 135°C (275°F) 3 min.	50 each/box (4 boxes/case)
1243A	3M™ Attest™ Steam Chemical Integrator (AAMI Type 5)¹	Dimensions: 5,1cm x 1,9cm (2 in x 3/4 in)	500 each/bag (2 bags/case)
1243B			100 each/bag (10 bags/case)
1355-18mm	3M™ Comply™ Lead Free Steam Indicator Tapes for Disposable Wraps (AAMI Type 1)¹	Dimensions: 18mm x 55m (0.70 in x 60 yds)	28 rolls/case
1355-24mm			20 rolls/case

Request your starter kit and save. It's everything you need to get started today.

490MKIT	3M™ Attest™ Super Rapid Mini Auto-reader Starter Kit	
<i>Kit includes:</i>	3M™ Attest™ Steam Chemical Integrator 1243B	100 each/bag (10 bags/case)
	3M™ Comply™ Lead Free Steam Indicator Tape for Disposable Wraps 1355-24mm	10 rolls/case
	3M™ Attest™ Super Rapid Readout Biological Indicator (Brown cap) 1492V*	50 each/box (2 boxes/case)
	3M™ Attest™ Mini Auto-reader 490M	1 unit/box



Contact your 3M authorized distributor representative to get started. Learn more at go.3M.com/490Mdentel.



Medical Solutions Division
3M Health Care
2510 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-228-3957
Web 3M.com/medical

*Refer to the indications for use provided in the biological indicator Instructions for Use. The indications for use should align with the sterilization cycle to be monitored.

¹ As defined by ANSI/AAMI/ISO 11140-1:2014.

© 2020 3M. All rights reserved. 3M, Attest and Comply are marks and/or registered marks of 3M. Unauthorized use prohibited. Please recycle. Printed in USA. EN_GL_70-2011-7947-3