



# Bair Hugger® Therapy Temperature Management Unit

Model 775



**Quiet. Adjustable. Safe. Effective.**  
It's everything you want from your forced-air warming unit...and more.



# Quiet, adjustable warming options for all Bair Hugger® blanket styles



Arizant Healthcare introduced forced-air warming to the world in 1987 with the launch of Bair Hugger therapy. Since that time, our Bair Hugger warming units and blankets have safely and effectively warmed more than 135 million patients worldwide.

Our latest advancement in forced-air warming – **the Bair Hugger Model 775 temperature management unit** – offers two airflow settings, hose-end temp sensing, precise temperature delivery and ease of use from a single, quiet warming unit. It's also compatible with all Bair Hugger blanket styles – pediatric to adult and Bair Paws® gowns.

Four temperature settings provide you the ability to deliver effective forced-air patient warming and cooling therapy. The ambient setting can complement your current patient cooling strategies.

Three temperature sensors and adjustable airflow provide safe, quick temperature response.

The Model 775 can track and report information on forced-air warming duration, over-temperature conditions and fault codes, all of which is accessible via the front panel – there's no need to open the unit.

Accessories for the Model 775 warming unit allow the unit to be mounted on a hard surface, an IV pole, bed rail or rolling cart.



©2010 Arizant Healthcare Inc. All rights reserved.  
10393 West 70th Street, Eden Prairie, MN 55344 USA  
tel 952-947-1200 | 800-733-7775  
fax 952-947-1400 | 800-775-0002  
web arizanthealthcare.com | bairhugger.com

Bair Hugger, Bair Paws, Arizant, Arizant Healthcare and the Bair Hugger and Arizant logos are trademarks of Arizant Healthcare Inc. registered in the U.S. Patent & Trademark Office.

603032C 4/10

## Specifications

### Dimensions

13'h x 14'w x 13'd (33h x 35w x 33d cm)

### Weight

16 lb (7.3 kg)

### Operating Temperatures

High:  $43^{\circ} \pm 1.5^{\circ}\text{C}$  ( $109.4^{\circ} \pm 2.7^{\circ}\text{F}$ )

Med:  $38^{\circ} \pm 1.5^{\circ}\text{C}$  ( $100.4^{\circ} \pm 2.7^{\circ}\text{F}$ )

Low:  $32^{\circ} \pm 1.5^{\circ}\text{C}$  ( $89.6^{\circ} \pm 2.7^{\circ}\text{F}$ )

Ambient: (Room Temperature)

### Fan Speeds and Relative Noise

High Fan Setting: 37 CFM (17.5 L/s) 53 dBA

Low Fan Setting: 32 CFM (15.1 L/s) 48 dBA

### Leakage Current

Meets UL 60601-1 and IEC 60601-1 requirements.

CAN/CSA-C22.2, No. 601.1

### Safety Systems

Audible and visual alarms;

Over-temperature  $\leq 56^{\circ}\text{C}$ ,  $53^{\circ}\text{C} \pm 3^{\circ}\text{C}$  typical

### Filter

High-efficiency 0.2  $\mu\text{m}$  filter

### Device Ratings

110-120 VAC, 50/60 Hz, 11.7 Amperes

220-240 VAC, 50/60 Hz, 7.2 Amperes

100 VAC, 50/60 Hz, 14.5 Amperes

### Compatibility

Use only with Bair Hugger blankets and Bair Paws gowns

## Ordering Information

Bair Hugger Model 775 Temperature Management Unit

**To order, please contact your Arizant Healthcare Inc. representative or call 1-800-733-7775.**