



**iM20**  
Patient Monitor

# iM20 Patient Monitor



## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Edan Instruments, Inc. | 3/F-B, Nanshan Medical Equipments Park,  
1019# Nanhai Rd., Shenzhen | 518067 P.R. China  
+86.755.26898326 | [www.edan.com.cn](http://www.edan.com.cn) | [info@edan.com.cn](mailto:info@edan.com.cn)

© 2015 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission.



ENG-PM-iM20-V3.0-20151218



# iM20

## Patient Monitor



Serving as both multi-parameter module of elite V series monitors and transport monitor, iM20 provides a comprehensive and seamless monitoring solution both inside and outside hospital, connecting ambulatory cares to intensive cares.

### > Parameter Configurations

- iSEAP™ ECG algorithm optimized for arrhythmia detection, pacemaker detection, and HR measurement
- SEMIP® 12-lead diagnosis algorithm with 208 kinds of ECG findings in optional 12-lead ECG configuration
- iMAT™ SpO<sub>2</sub> algorithm with outstanding motion resistance and low perfusion resistance performance
- Perfusion Index (PI) measurement as a blood perfusion indicator
- Optional Nellcor SpO<sub>2</sub> with SatSeconds™ technology
- iCUFFS™ NIBP algorithm optimized for cardiac patients, hypertensive patients, and neonatal patients
- Optional SunTech NIBP with SunTech Advantage TMT technology
- Optional 2-channel invasive blood pressure with waveform overlapping function
- Optional EFM module with Respironics LoFlo™ technology inside
- 2-channel temperature monitoring



### > Built-in Temporary Memory

150 hours	1200 sets	200 sets	200 sets	50 sets
Trend Review	NIBP Measurement Review	Alarm Review	Arrhythmia Review	12-lead Diagnosis Review

- Portable Design
- Full Touch Operation
- IP44 Waterproof Dustproof
- 1.2 meters Shock Resistance
- 4 hours Rechargeable Li-ion Battery
- No-fan Design
- 256M Built-in Memory

### > Ambulatory Designs



Designed as a robust transport monitor, iM20 is a solid companion in the harshest conditions.

- With IP44 water and dust protection, iM20 is outstanding for field medical practice, from regular emergency cares in rainy weather, to disaster controls with dusty environments.
- Designed with 1.2 meter shock resistance, unintentional falls and impacts shall not interrupt the monitoring, nor does the bumpy ride on rough road or hilly areas.
- With the ability to resist extreme temperature and altitudes, iM20 can be used in cold mountain areas as well as hot tropical places.

