


VE-100/VE-300

Veterinary Electrocardiograph




EDAN

A world of potential

VE-100/VE-300 Veterinary Electrocardiograph

We respect lives, and treat them equally. Our Three-Channel and Single-Channel ECG machines for vets are carefully designed without compromising in quality and accuracy, ensuring your confidence in dealing with your lovely pets.



VE-100 Single-Channel ECG

- Compact and portable design
- Vivid waveform display
- Built-in rechargeable Li-ion battery, AC/DC power supply
- Auto/Manual work mode selectable
- Lead-off detection and display
- Complete digital filters, resisting baseline drift, AC and EMG interference
- Baseline automatically adjustable to optimize printing position
- High-resolution thermal printing

VE-300 Three-Channel ECG

- Compact and portable design
- Vivid waveform display
- Built-in rechargeable Li-ion battery, AC/DC power supply
- Auto/Manual work mode selectable
- Lead-off detection and display
- Complete digital filters, resisting baseline drift, AC and EMG interference
- Baseline automatically adjustable to optimize printing position
- High-resolution thermal printing
- Communication interface: Ethernet/RS232/USB
- PC-based ECG data management software (optional)

