

## **EMERGENCY - SQ SQUARE TROLLEY** **Q240 W.(5)CR**

**DIMENSIONS** 480 x 480 mm

**ADEQUATE HEIGHT, ERGONOMICS** 980 mm (work surface – floor).

### **WORK SURFACE**

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

### **MONOBLOC STRUCTURE**

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

### **SQ SQUARE DRAWERS**

- Polypropylene, PP. 2 colours - B Blue, W White.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 10 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

### **EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION**

- Elevated push handle. Ergonomics, height 1020 mm.
- Double wheels Ø 100 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at push handle side.
- Wheels Ø 125 mm ref. AD12 (optional).



## **Q240 W.(5)CR**

**SQ065 W.AST.12.14.30.34.36.40.(5)CR**

### **ACCESSORIES**

**AST** ELEVATED PUSH HANDLE

**12** MANUAL LOCKING

**14** NUMBERED DISPOSABLES SEALS - 100 units

**30** FOLD DOWN TRAY

**34** CARDIAC MASSAGE BOARD

**36** GAS BOTTLE STAND

**40** SERUM HOLDER

**(5)CR** RED CARDS