



CE 0123

MINIMAX is a device for aerosoltherapy. Designed for home-care with compact size and carrying bag for easy transport of unit, accessories and medications.

STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
- Adult Masks with elastic strings
- Child Masks with elastic strings
- Mouthpiece
- Nosepiece (non invasive)
- Air Tube
- Air Filter
- Carrying Bag
- Detachable power cord

TECHNICAL FEATURES

Compressor	Oilless and maintenance-free piston pump
Power Feeding	230V-50Hz - other configurations on request
Power consumption	170 VA
Max Pressure	2.5 bar 36 psi 250 kPa
Operating Pressure	1.1 bar 16 psi 110 kPa
Max Air Flow	14 l/min
Operating Air Flow	4,8 l/min
Neb-rate with 4 ml solution	0,35 ml/min
Mass Median Aerodynamic Diameter	2,44 μ m (MMAD) according to EN-13544-1
Respirable Fraction (< 5 μ m)	> 80%
Duty Cycle	Ton: 20min Toff: 40min
Noise Level (EN 13544-1 Standard)	Approx. 60 dBA
Weight	1,45 kg
Size (Main Unit)	cm 13 x 10 x 15
Years of Warranty	3
Shipping carton	6
Place of Manufacturing	Italy



Italian
Medical
Touch
MADE IN ITALY

CM005.REV1 - EDITION NOVEMBER 2014

CODE

RE 300250 ————— MINIMAX