


FEATURES OF JBP NANONEEDLE



Ultra-Thin Walled


Extra Fine

Micron Polishing Technology
(Micron Sharpness)



- Less Pain & Swelling
- Little Bleeding
- Invisible Injection Mark
- Injectable For High Viscid Solution

Comparison of Inner and Outer Diameters between JBP Nanoneedle and Conventional Needle



Gauge	JBP Nanoneedle		Conventional Needle	
	Outer Diameter	Inner Diameter	Outer Diameter	Inner Diameter
34G	0.20 mm	0.10 mm	-	-
33G	0.26 mm	0.16 mm	0.26 mm	0.10 mm
30G	0.30 mm	0.20 mm	0.30 mm	0.15 mm

Comparison of Injection Marks between JBP Nanoneedle and Conventional Needle

