


Opal Examination Bed



Item Code: 003036

Width 闊度: 610mm, Length 長度: 1910mm

Height adjustment 高度: 495mm - 995mm

Back rest tilt 背板斜度: 75°

Foot rest tilt 腳板斜度: 50°

Fire Retardant 防火物料

Lifting capacity 承重: 260kg

Wheels 車輪: 125mm with brakes

Colour 顏色: Blue

Optional 另可加配
Extra Width- 710mm
Foot switch
Cardio cut

Set includes

Paper roll holder (AMC PRH)