

3M™ Surgical Drapes, Gowns and Custom Procedure Trays Range Guide



- Single Use Drapes
- Surgical Drape Packs
- Custom Procedure Trays
- Gowns
- Bespoke Soft Packs
- Incise Drapes



The 3M™ Surgical Drape range...

From the leader in Medical Adhesive and Drape Materials

Surgeons all over the world rely on 3M for many of their draping needs. 3M has made use of innovative technologies to create a highly advanced range of single-use surgical drapes, for the highest level of performance and reliability.

Cost Savings

The 3M surgical drape range consists of both specialised packs and individual drapes for cost effective and flexible draping solutions.

Customer Support

3M provides a high-quality and dependable level of service to thousands of health care workers each day. This means offering you the very best in customer service; a fully dedicated team of experienced theatre nurses will assist you in meeting the fast changing demands of today's surgical procedures.

Unique Drape Exchange Scheme

We recognise that lists and procedures are subject to last minute changes. Therefore your 3M representative will be happy to exchange any uncontaminated opened and unused 3M surgical drape or surgical gown, free of charge, for a new sterile drape or gown.



...the role of Surgical Drapes

Surgical Drapes are designed to establish and maintain a sterile field. It is very important to select a surgical drape that provides reliable performance. To achieve the best surgical outcomes, look for surgical drapes with these qualities:

3M™ Steri-Drape™ 9000 Surgical Drapes

Steri-Drape 9000 material consists of a laminate of polypropylene fibres and a polyethylene film. These materials are by-products of the oil industry. There is no cellulose in these drapes, which means that no trees are used in their manufacture. The products do not present a significant risk to the environment when incinerated after use and are broken down to carbon dioxide and water.

Impervious to Strikethrough

A drape that allows fluid transfer, and thus bacteria transfer, compromises the sterile field, which increases the risk of wound infection. That is why a drape that is impervious to liquid strikethrough is so important.

Fluid Control Features

3M surgical drapes effectively control and contain fluid. This not only reduces potential exposure to infectious body fluids, but keeps the patient dry which can help reduce turnaround time between procedures. Many 3M surgical drapes have fluid collection pouches, and the 3M™ Steri-Drape™ 9000 Range of drapes is absorbent.

Reliable Adhesion



3M is a leader in adhesive technology, and has applied this expertise to provide surgical drapes with consistently reliable adhesion to the skin or wound edge. Even during heavy retraction and irrigation, 3M incise drapes adhere securely to wound edges, providing a barrier to bacteria on the skin.

Superior Strength and Comfort

The material of the Steri-Drape 9000 range of surgical drapes is strong enough to withstand the most arduous procedures, and even the heaviest surgical instruments can be placed on them without concern. They are also soft and pliable for both patient and staff.

3M Range

3M surgical drapes are available individually and in specialised packs to cover most surgical procedures. The contents of each draping pack are tailored to match the requirements of each and every procedure. Fully customised packs are also available, to save time and minimise waste.



3M Service and Support

Clinical Nurse Advice

Our people; your support team

3M provides a high quality and dependable level of service to thousands of health care workers each day. This means offering all the support possible to help meet the fast changing demands of today's surgical procedures.

3M's commitment to providing the very best in customer service is reflected in the fact that we have a fully dedicated team of experienced theatre nurses in addition to your account manager, who ensure that we provide on-going support to our customers. The team can support you with in-house training and can also deliver a number of training packages.

Behind your immediate 3M account support team, there exists a full community of technical and quality assurance employees. Our UK based, experienced team of Customer Services staff are very experienced in the Health Care market. We also have a Technical Services Team who are available 24 hours a day to cover technical emergencies out of hours.

3M Environmental Care

3M has an approach to the environment based on prevention at source.

The company believes that it is better to prevent pollution rather than remove it after it has been created. This concept means applying more intelligence to produce high-value products with fewer resources and fewer environmental impacts.

The 3M environmental policy involves:

- Preventing pollution at source, whenever possible.
- Developing products that have a minimum effect on the environment.
- Conserving natural resources through reclamation and other appropriate methods.
- Meeting and sustaining government regulations.
- Participating in government-led environmental initiatives, wherever possible.

With strong environmental management systems in place, 3M has a long history of continuous improvements to minimise the impact of our manufacturing processes and products on the global environment. Over the last 6 years our global environmental progress has included a 26% improvement in energy efficiency, 25% reduction in waste, 58% reduction in volatile organic air emissions and the successful completion of 895 'Pollution Prevention Pays' projects.

Life Cycle Management (LCM) is an integral part of our environmental activities. Instead of only focusing on the manufacturing process to control environmental, health, safety and energy effects, the spotlight is now on products throughout their entire life cycle, from manufacturing through to customer use and disposal. 3M has adopted a Life Cycle Management Policy requiring all new and existing products to undergo a LCM review, which must be completed by 2010.

In the UK and globally, 3M has been recognised for its commitment to the environment through a number of prestigious awards, including both a Stratospheric Ozone Protection Award and a 2006 Energy Star Award for Sustained Excellence in Energy Management from the US Environmental Protection Agency.

Selection Guide

3M™ Custom Procedure Trays	10
3M™ Incise Drapes	12-16
3M™ Ioban™ 2 Antimicrobial Incise Drapes	13
3M™ Ioban™ EZ Antimicrobial Incise Drapes	14
3M™ Steri-Drape™ Incise Drapes	15
3M™ Steri-Drape™ 2 Incise Drapes	16
3M™ General Surgery Drapes	18-26
3M™ Orthopaedic Drapes	28-36
3M™ Cardiovascular Drapes/Cardiac Catheterisation Drapes	38-40
3M™ Obstetrics and Gynaecology Drapes	42-48
3M™ Urology Drapes	50-54
3M™ Ophthalmic Drapes	56-60
3M™ Head and Neck Drapes	62-64
3M™ Neurosurgery Drapes	66-67
3M™ Accessories and other Drapes	70-82
3M™ Disposable Surgical Gowns Standard and Reinforced	84-86
3M™ High Performance Surgical Gowns – Standard and Reinforced	84-85
3M™ Disposable Surgical Gowns – Standard and Reinforced	84-86
Cross reference by product number	88-89

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

Cardiovascular Drapes/
Cardiac Catheterisation Drapes

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

Disposable Surgical Gowns
– Standard and Reinforced

Cross Reference by Product Number

3M™ Custom Procedure Trays

3M™ Custom Procedure Trays

– *Reliable; Flexible; Cost effective*

Whatever your requirements, 3M Health Care can design a procedure tray that meets your needs. We can access over 20,000 components from other leading health care manufacturers, and you have our assurance that every pack will still conform to our high standards of quality, reliability and performance. Please contact our Customer Service Team on 01509 613151 for further details of how our customised service can help you.

In addition to customised procedure trays, standard drape packs are also available. Please see the relevant procedure sections for details.



Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

**Cardiovascular Drapes/
Cardiac Catheterisation Drapes**

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

**Disposable Surgical Gowns
– Standard and Reinforced**

Cross Reference by Product Number

3M™ Incise Drapes

3M™ Ioban™ 2 Antimicrobial Incise Drapes

- Contain iodophor which is slowly released on to the skin surface giving continuous and controlled antimicrobial activity throughout the surgical procedure.
- Have a high moisture vapour transmission rate which allows the skin to “breathe” ensuring improved adhesion throughout the longest procedures.
- Are highly conformable to body contours. The unique low memory stretch allows limb mobilisation or heavy retraction with reduced tension to the skin.
- Have superb adhesion right up to the wound edge.



3M™ Ioban™ EZ Antimicrobial Incise Drapes

- In addition, the Ioban EZ incise drape is also available which includes the “easy-delivery” drape application feature (except 6635), including a removable cardboard strip for greater stability.



3M™ Steri-Drape™ 2 Incise Drapes

- Are made of a breathable film which provides superior adhesion.
- Have low memory stretch.
- Conform to body contours and stretch to allow limb manipulation.
- Reduce or eliminate the need for towel clips by stabilising surgical drape sheets.




3M™ Steri-Drape™ Incise Drapes

- Are the original incise drapes which are totally impermeable, have a pressure sensitive adhesive and provide a sterile operating surface.
- Reduce or eliminate the need for towel clips by stabilising surgical drape sheets.



3M™ Ioban™ 2 Antimicrobial Incise Drapes



6635 Antimicrobial Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
15cm x 20cm	10cm x 20cm	10	4



6640 Antimicrobial Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
44cm x 35cm	34cm x 35cm	10	4




6648 Antimicrobial Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
66cm x 60cm	56cm x 60cm	10	4



6650 Antimicrobial Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
66cm x 45cm	56cm x 45cm	10	4




6651 Antimicrobial Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
66cm x 85cm	56cm x 85cm	10	4


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Ioban™ EZ Antimicrobial Incise Drapes




6640EZ Antimicrobial Incise Drape

Incise Area	Drapes per Box	Boxes per Case
35cm x 35cm	10	4




6648EZ Antimicrobial Incise Drape

Incise Area	Drapes per Box	Boxes per Case
60cm x 60cm	10	4



6650EZ Antimicrobial Incise Drape

Incise Area	Drapes per Box	Boxes per Case
60cm x 45cm	10	4




6651EZ Antimicrobial Incise Drape

Incise Area	Drapes per Box	Boxes per Case
60cm x 85cm	10	4


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Steri-Drape™ Incise Drapes




1035 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
15cm x 20cm	10cm x 20cm	10	4




1037 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
29cm x 25cm	19cm x 25cm	10	4




1040 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
60cm x 35cm	35cm x 35cm	10	4



1050 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
82cm x 45cm	60cm x 45cm	10	4




1051 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
82cm x 85cm	60cm x 85cm	10	4


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Steri-Drape™ 2 Incise Drapes




2035 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
20cm x 20cm	10cm x 20cm	10	4




2037 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
38cm x 25cm	28cm x 25cm	10	4




2040 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
38cm x 41cm	28cm x 41cm	10	4




2045 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
60cm x 45cm	50cm x 45cm	10	4



2050 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
90cm x 60cm	60cm x 60cm	10	4



2051 Incise Drape

Overall Size	Incise Area	Drapes per Box	Boxes per Case
90cm x 85cm	60cm x 85cm	10	4

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

**Cardiovascular Drapes/
Cardiac Catheterisation Drapes**

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

**Disposable Surgical Gowns
– Standard and Reinforced**

Cross Reference by Product Number

General Surgery Drapes

As general surgery procedures vary so much the drapes used need to be adaptable and cater for the challenges likely to be met, such as large volumes of blood and irrigation fluids, heavy retraction during major abdominal procedures and procedures that may last for several hours.

Impervious to Strikethrough

- A drape that allows fluid transfer, and thus bacteria transfer, compromises the sterile field, which increases the risk of wound infection. That is why a drape that is impervious to liquid strikethrough is so important.

Fluid Control Features

- 3M surgical drapes effectively control and contain fluid. This not only reduces potential exposure to infectious body fluids, but keeps the patient dry which can help reduce turnaround time between procedures. Many 3M surgical drapes have fluid collection pouches, and the 3M™ Steri-Drape™ 9000 Range of drapes is absorbent.

Superior Strength and Comfort

- The material of the Steri-Drape 9000 range of surgical drapes is strong enough to withstand the most arduous procedures, and even the heaviest surgical instruments can be placed on them without concern. They are also soft and pliable for both patient and staff.

3M™ Absorbent Impervious Fabric is:


- Low-linting, non-woven, absorbent fabric laminated to a plastic film
- Prevents fluids from running off the drape surface
- Absorbent and impervious throughout the entire drape
- Impervious to strikethrough
- Fluid control features
- Superior strength and comfort

3M™ Ioban™ 2 Antimicrobial Incise Drapes

- Ioban 2 antimicrobial incise drapes contain iodophor which is slowly released onto the skin surface giving continuous and controlled antimicrobial activity throughout the surgical procedure. An Ioban incise area is featured in some speciality drapes. Ioban incise drapes are also available as separate items.




3M™ General Surgery Drapes




9000 Basic Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Adhesive Towel Drapes	2	100cm x 75cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		



9000L Basic Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	5	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Adhesive Towel Drapes	2	100cm x 75cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (XL)	1	300cm x 228cm		
Instrument Table Cover	1	200cm x 150cm		




9007 ENT Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	2	10cm x 55cm	8	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Drape Sheet	1	75cm x 90cm		
Adhesive Split Sheet	1	150cm x 250cm		
Adhesive Towel Drape	1	100cm x 75cm		
Instrument Table Cover	1	150cm x 150cm		


KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	

3M™ General Surgery Drapes




9008 Head and Neck Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Headwrap	1	114cm x 90cm		
Adhesive Split Sheet	1	150cm x 250cm		
Instrument Table Cover	1	200cm x 150cm		



9010 Basic Surgical Pack without Mayo

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Adhesive Towel Drapes	2	100cm x 85cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		




9018 Surgical Split Sheet Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	2	10cm x 55cm	6	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Adhesive Towel Drape	1	100cm x 75cm		
Adhesive Split Sheet	1	228cm x 260cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ General Surgery Drapes




9001 Paediatric Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Adhesive Towel Drapes	2	75cm x 75cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		
Narrow Adhesive Strip				




9040 Leggings (pair)

Overall Size	Drapes per Box	Boxes per Case
75cm x 114cm	12	4



9050 Adhesive Aperture Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
228cm x 300cm	12cm x 16cm	10	2




9051 Thyroid Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
183cm x 315cm	12cm x 12cm	15	2


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ General Surgery Drapes




9015 Hand Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Hand Drape	1	360cm x 200cm	7	2
Hand Towels	4	-		
Mayo Stand Cover	1	80cm x 145cm		
Instrument Table Cover	1	200cm x 150cm		




9062 Laparotomy Drape (Transverse)

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
254cm x 285cm	30cm x 10cm	12	2



9063 Laparotomy Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
254cm x 285cm	10cm x 30cm	10	2




1092 Minor Procedure Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
56cm x 64cm	6.3cm	25	4


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ General Surgery Drapes



1115 Tube Organiser


Overall Size	Units per Box	Boxes per Case
14cm x 4cm	100	4



1016 Irrigation Pouch

Overall Size	Units per Box	Boxes per Case
50cm x 55cm	10	4


- Adhesive Strip
- 1 Exit Port



1084 Under Buttocks Drape with Pouch

Overall Size	Drapes per Box	Boxes per Case
86cm x 108cm	30	4


- Fluid Collection Pouch



2057 Incise Pouch

Overall Size	3M™ Steri-Drape™ 2 Incise Area	Drapes per Box	Boxes per Case
87cm x 74cm	30cm x 30cm	10	4

- 360-Degree Fluid Collection Pouch
- Tube Organisers




2058 Aperture Pouch

Overall Size	Fenestrated Incise Area	Drapes per Box	Boxes per Case
87cm x 74cm	30cm x 30cm	10	4

- 360-Degree Fluid Collection Pouch
- Tube Organisers

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	


3M™ General Surgery Drapes



6657 Pouch with 3M™ Ioban™ 2 Incise Film

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
87cm x 74cm	30cm x 30cm	10	4


- 360-Degree Fluid Collection Pouch
- Tube Organisers



6658 Large Pouch with 3M™ Ioban™ 2 Incise Film

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
74cm x 74cm	33cm x 43cm	5	4


- Fluid Collection Pouch



6659 Extra Large Pouch with 3M™ Ioban™ 2 Incise Film

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
74cm x 87cm	42cm x 52cm	5	4


- 2 Exit Ports



1003 Isolation Bag

Overall Size	Units per Box	Boxes per Case
49cm x 49cm	10	4

- Drawstring Closure
- Clear Plastic




Wound Edge Protector

	Overall Size	Ring Diameter	Packs per Box	Boxes per Case
1073	90cm x 90cm	12.5cm	10	4
1074	90cm x 90cm	17.5cm	10	4
1075	90cm x 90cm	22.5cm	10	4
1076	90cm x 90cm	27.5cm	10	4
1077	90cm x 90cm	9.5cm	10	4


- 4 Adhesive Patches
- 1 Flexible Ring

3M™ General Surgery Drapes




9052 Hand Drape

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
360cm x 200cm	3cm	9	2



9054 Foot / Hand Drape

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
200cm x 300cm	3cm	12	2



2065 3M™ Abdomino-Perineal Drape-Lithotomy Position

Overall Size	3M™ Steri-Drape™ 2 Incise Area	Pads per Box	Boxes per Case
195cm x 271cm	41cm x 45cm	8	2


- Tube Organisers
- Perineal Adhesive Aperture: 12cm x 15.5cm

2065P Abdomino-Perineal Draping Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Abdomino-Perineal Drape	1	195cm x 271cm	5	2
Instrument Table Cover	1	100cm x 150cm		
Mayo Stand Cover	1	80cm x 145cm		
Towel Drape	1	75cm x 90cm		
Hand Towels	4	-		

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ General Surgery Drapes




6665 3M™ Ioban™ 2 Abdomino-Perineal Drape - Lithotomy Position

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
200cm x 270cm	35cm x 46cm	10	1

- Tube Organisers
- Perineal Adhesive Aperture: 12cm x 15.5cm

6665P001 Abdomino-Perineal Draping Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Abdomino-Perineal Drape	1	200cm x 270cm	10	1
Instrument Table Cover	1	100cm x 150cm		
Mayo Stand Cover	1	80cm x 145cm		
Towel Drape	1	75cm x 90cm		
Hand Towels	4	-		




1032 Laparoscopic Drape - Lithotomy Position

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
241cm x 228cm	27cm x 30cm	8	2

- 4 Tube Organisers (ref. 1115)
- Perineal Adhesive Aperture: 12cm x 15.5cm
- Available as a single box

1032P Laparoscopic Drape Pack / Lithotomy Position

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Laparoscopy Drape	1	241cm x 228cm	5	2
Instrument Table Cover		200cm x 150cm		
Hand Towels	4	-		
Mayo Stand Cover	1	80cm x 145cm		



1031 Laparoscopic Drape - Supine Position

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
285cm x 182cm	27cm x 30cm	10	2

- 4 Tube Organisers (ref. 1115)
- Available as a single box

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

Cardiovascular Drapes/
Cardiac Catheterisation Drapes

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

Disposable Surgical Gowns
– Standard and Reinforced

Cross Reference by Product Number

3M™ Orthopaedic Drapes

Orthopaedic surgery requires a draping system that maintains a sterile field throughout the procedure and reliably adheres to the skin to reduce the likelihood of drapes moving over long periods of time. Heavy manipulation often occurs during orthopaedic procedures and drapes need to be able to withstand this.

Superior Strength and Comfort

- The material of the 3M™ Steri-Drape™ 9000 Range of surgical drapes is strong enough to withstand the most arduous procedures, and even the heaviest surgical instruments can be placed on them without concern. They are also soft and pliable for both patient and staff.

Impervious to Strikethrough

- A drape that allows fluid transfer, and thus bacteria transfer, compromises the sterile field, which increases the risk of wound infection. That is why a drape that is impervious to liquid strikethrough is so important.

Reliable Adhesion

- 3M is a leader in adhesive technology, and has applied this expertise to provide surgical drapes with consistently reliable adhesion to the skin or wound edge. Even during heavy retraction and irrigation, 3M incise drapes adhere securely to wound edges, providing a barrier to bacteria on the skin.

Fluid Control Features


- 3M surgical drapes effectively control and contain fluid. This not only reduces potential exposure to infectious body fluids, but keeps the patient dry which can help reduce turnaround time between procedures. Many 3M surgical drapes have fluid collection pouches, and the Steri-Drape 9000 range of drapes is absorbent.

3M™ Ioban™ 2 Antimicrobial Incise Drapes

- Ioban 2 antimicrobial incise drapes contain iodophor which is slowly released onto the skin surface giving continuous and controlled antimicrobial activity throughout the surgical procedure. An Ioban incise area is featured in some speciality drapes. Ioban incise drapes are also available as separate items.




3M™ Orthopaedic Drapes




9002 Orthopaedic Surgical Pack without Mayo

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	3	10cm x 55cm	5	2
Hand Towels	4	–		
Legging	1	35cm x 100cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Split Sheet	1	228cm x 260cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		



9004 Extremity Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	2	10cm x 55cm	5	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Drape Sheet	1	114cm x 150cm		
Adhesive Towel Drape	1	100cm x 75cm		
Extremity Drape	1	228cm x 300cm		
Instrument Table Cover	1	200cm x 150cm		




9000 Basic Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Adhesive Towel Drapes	2	100cm x 75cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Orthopaedic Drapes




9012 Orthopaedic Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	3	10cm x 55cm	5	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Legging	1	35cm x 100cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Split Sheet	1	228cm x 260cm		
Adhesive Drape Sheet (L)	1	150cm x 250cm		
Instrument Table Cover	1	200cm x 150cm		



9012L Orthopaedic Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	3	10cm x 55cm	4	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Legging	1	35cm x 100cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Split Sheet	1	228cm x 260cm		
Adhesive Drape Sheet (XL)	1	300cm x 228cm		
Instrument Table Cover	1	200cm x 150cm		




9014 Extremity Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	2	10cm x 55cm	5	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Drape Sheet	1	114cm x 150cm		
Legging	1	35cm x 60cm		
Extremity Drape	1	228cm x 300cm		
Instrument Table Cover	1	200cm x 150cm		


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Orthopaedic Drapes




9010 Basic Surgical Pack without Mayo

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Adhesive Towel Drapes	2	100cm x 85cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		




Adhesive Split Sheets

	Overall Size	Drapes per Box	Boxes per Case
Split Size 7 x 70cm 9043	150cm x 183cm	30	2
9044	150cm x 250cm	25	2
Split Size 20 x 100cm 9045	228cm x 260cm	15	2



9052 Hand Drape

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
360cm x 200cm	3cm	9	2




9015 Hand Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Hand Drape	1	360cm x 200cm	7	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Instrument Table Cover	1	200cm x 150cm		


KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	

3M™ Orthopaedic Drapes




9053 Extremity Drape

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
228cm x 300cm	6cm x 8cm	10	2



9054 Foot / Hand Drape

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
200cm x 300cm	3cm	12	2




9055 One-Piece Hip Drape with Side Pouches

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
300cm x 353cm	43cm x 97cm	8	2

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Orthopaedic Drapes




9063 Laparotomy Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
254cm x 285cm	10cm x 30cm	10	2



Leggings


	Overall Size	Drapes per Box	Boxes per Case
9041	35cm x 60cm	50	4
9042	35cm x 100cm	35	4



1014 Patient Isolation Drape

Overall Size	3M™ Steri-Drape™ Incise Area	Drapes per Box	Boxes per Case
320cm x 213cm	48cm x 25cm	5	4

- Adhesive Strip Along Top
- Clear Plastic



1017 Isolation Drape with Incise Film and Pouch

Overall Size	3M™ Steri-Drape™ Incise Area	Drapes per Box	Boxes per Case
320cm x 213cm	48cm x 25cm	5	4


- Fluid Collection Pouch
- 1 Exit Port
- 2 Instrument Pouches
- 2 Tube and Cord Organisers
- Adhesive Strip Along Top
- Clear Plastic

1017P001 Patient Isolation Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Patient Isolation Drape	1	320cm x 213cm	12	1
Instrument Table Cover	1	200cm x 150cm		
Mayo Stand Cover	1	80cm x 145cm		
Hand Towels	4	-		

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Orthopaedic Drapes




130cm

1015 U-Drape

Overall Size	Drapes per Box	Boxes per Case
120cm x 130cm	10	4

- U-Slot Aperture with Adhesive
- Clear Plastic




180cm

1019 Long U-Drape

Overall Size	Drapes per Box	Boxes per Case
120cm x 180cm	10	4

- U-Slot Aperture with Adhesive
- Adhesive Patch
- Clear Plastic



213cm


6617 Isolation Drape with 3M™ Ioban™ 2 Incise Film and Pouch

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
320cm x 213cm	48cm x 25cm	5	4

- Fluid Collection Pouch
- 1 Exit Port
- 2 Instrument Pouches
- 2 Tube and Cord Organisers
- Adhesive Strip Along Top
- Clear Plastic

6617P001 Patient Isolation Pack with 3M™ Ioban™ 2 Incise Film and Pouch

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Patient Isolation Drape	1	320cm x 213cm	14	1
Instrument Table Cover	1	200cm x 150cm		
Mayo Stand Cover	1	80cm x 145cm		
Hand Towels	4	-		



252cm


6619 Large Isolation Drape with 3M™ Ioban™ 2 Incise Film and Pouch

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
328cm x 252cm	70cm x 32cm	5	1

- Fluid Collection Pouch
- 1 Exit Port
- 2 Instrument Pouches
- 2 Tube and Cord Organisers
- Adhesive Strip Along Top
- Clear Plastic

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area


3M™ Orthopaedic Drapes



1094 Arthroscopy Sheet with Pouch

Overall Size	Drapes per Box	Boxes per Case
226cm x 320cm	20	1

- Fluid Collection Pouch
- 1 Exit Port
- Tube and Cord Organiser




1194 Arthroscopy Sheet with Pouch

Overall Size	Drapes per Box	Boxes per Case
226cm x 320cm	5	1

- Fluid Collection Pouch
- 2 Exit Ports
- Tube and Cord Organiser

1195 Arthroscopy Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Arthroscopy Sheet with Pouch	1	226cm x 320cm	5	1
Long U-drape	1	120cm x 180cm		
Impervious Stockinette	1	23cm x 71cm		
3M™ Coban™ Self Adherent Bandage	1	3 metres x 10cm		
Mayo Stand Cover	1	55cm x 135cm		
Back Table Cover	1	120cm x 225cm		




9196 Shoulder Split Sheet with Pouch

Overall Size	Drapes per Box	Boxes per Case
228cm x 260cm	10	2

- Fluid Collection Pouch
- 2 Exit Ports

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area


3M™ Orthopaedic Drapes



1007 Arthroscopy Fluid Collection Pouch

Overall Size	Drapes per Box	Boxes per Case
115cm x 58cm	10	4


- Adhesive Strips
- 1 Exit Port



1016 Irrigation Pouch


Overall Size	Drapes per Box	Boxes per Case
50cm x 55cm	10	4

- Adhesive Strip
- 1 Exit Port



1115 Tube Organiser


Overall Size	Drapes per Box	Boxes per Case
14cm x 4cm	100	4



1067 U-Pouch

Overall Size	Drapes per Box	Boxes per Case
87cm x 74cm	5	4

- U-Slot Aperture with Adhesive
- 2 Exit Ports



Stockinette

	Overall Size	Drapes per Box	Boxes per Case
9109	23cm x 71cm	20	1
9110	30cm x 122cm	20	1

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

**Cardiovascular Drapes/
Cardiac Catheterisation Drapes**

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

**Disposable Surgical Gowns
– Standard and Reinforced**

Cross Reference by Product Number

3M™ Cardiovascular Drapes / Cardiac Catheterisation Drapes

Most cardiothoracic procedures require drapes that can withstand heavy manipulation and reduce the likelihood of infection. There is often a large amount of fluid involved with these procedures and therefore drapes used should have fluid control features. Drapes that contain antimicrobial features can reduce the risk of surgical site infection which is very important in this type of surgery.



Superior Strength and Comfort

- The material of the 3M™ Steri-Drape™ 9000 Range of surgical drapes is strong enough to withstand the most arduous procedures, and even the heaviest surgical instruments can be placed on them without concern. They are also soft and pliable for both patient and staff.



Impervious to Strikethrough

- A drape that allows fluid transfer, and thus bacteria transfer, compromises the sterile field, which increases the risk of wound infection. That is why a drape that is impervious to liquid strikethrough is so important.


Fluid Control Features

- 3M surgical drapes effectively control and contain fluid. This not only reduces potential exposure to infectious body fluids, but keeps the patient dry which can help reduce turnaround time between procedures. Many 3M surgical drapes have fluid collection pouches, and the Steri-Drape 9000 range of drapes is absorbent.




3M™ Ioban™ 2 Antimicrobial Incise Drapes

- Ioban 2 antimicrobial incise drapes contain iodophor which is slowly released onto the skin surface giving continuous and controlled antimicrobial activity throughout the surgical procedure. An Ioban incise area is featured in some speciality drapes. Ioban incise drapes are also available as separate items.




3M™ Cardiovascular Drapes / Cardiac Catheterisation Drapes



6677 Cardiovascular Sheet with 3M™ Ioban™ 2 Incise Film


Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
274cm x 338cm	41cm x 92cm	6	1

- Dual Side Pockets




9046 Bilateral Split Sheet

Overall Size	Aperture Split	Drapes per Box	Boxes per Case
183cm x 200cm	15cm x 35cm	15	2



9007 ENT Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	2	10cm x 55cm	8	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Drape Sheet	1	75cm x 90cm		
Adhesive Split Sheet	1	150cm x 250cm		
Adhesive Towel Drape	1	100cm x 75cm		
Instrument Table Cover	1	150cm x 150cm		




1090 Femoral Angiography Drape

Overall Size	Oval Aperture Split	Drapes per Box	Boxes per Case
182cm x 228cm	6cm x 9cm	20	1

- 2 Adhesive Patches: 11.5cm x 14cm

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	


3M™ Cardiovascular Drapes / Cardiac Catheterisation Drapes



6681 3M™ Ioban™ 2 Cardiovascular Drape


Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
230cm x 230cm	See Illustration	8	1

- Impervious Fabric
- Plastic Isolation Screen



9018 Surgical Split Sheet Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	2	10cm x 55cm	6	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Adhesive Towel Drape	1	100cm x 75cm		
Adhesive Split Sheet	1	228cm x 260cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		



1191 Femoral Angiography Drape

Overall Size	Oval Aperture Split	Drapes per Box	Boxes per Case
228cm x 330cm	6cm x 9cm	20	1

- Highly Absorbent Area: 91.4cm x 122cm
- 2 Adhesive Patches: 11.5cm x 14cm
- Transparent Side and End Panels

For X-ray equipment covers see page 81.

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

Cardiovascular Drapes/
Cardiac Catheterisation Drapes

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

Disposable Surgical Gowns
– Standard and Reinforced

Cross Reference by Product Number

3M™ Obstetrics and Gynaecology Drapes

Some obstetric procedures such as caesarean section present challenges for a draping system. Large amounts of blood and amniotic fluid present mean a drape must have excellent fluid control features to ensure the patient remains dry throughout the procedure and reduces the likelihood of fluid spillages that could be a health and safety hazard to healthcare workers. Ease of application is imperative as these procedures are often emergencies.

Fluid Control Features

- 3M surgical drapes effectively control and contain fluid. This not only reduces potential exposure to infectious body fluids, but keeps the patient dry which can help reduce turnaround time between procedures. Many 3M surgical drapes have fluid collection pouches, and the 3M™ Steri-Drape™ 9000 Range of drapes is absorbent.

Superior Strength and Comfort

- The material of the Steri-Drape 9000 range of surgical drapes is strong enough to withstand the most arduous procedures, and even the heaviest surgical instruments can be placed on them without concern. They are also soft and pliable for both patient and staff.

Impervious to Strikethrough

- A drape that allows fluid transfer, and thus bacteria transfer, compromises the sterile field, which increases the risk of wound infection. That is why a drape that is impervious to liquid strikethrough is so important.


3M Range

- 3M surgical drapes are available individually and in specialised packs to cover most surgical procedures. The contents of each draping pack are tailored to match the requirements of each and every procedure.
- Fully customised packs are also available, to save time and minimise waste.

3M™ Ioban™ 2 Antimicrobial Incise Drapes


- Ioban 2 antimicrobial incise drapes contain iodophor which is slowly released onto the skin surface giving continuous and controlled antimicrobial activity throughout the surgical procedure. An Ioban incise area is featured in some speciality drapes. Ioban incise drapes are also available as separate items.






9000 Basic Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Adhesive Towel Drapes	2	100cm x 75cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		



9003 Cystoscopy Surgical Pack


Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Drape Sheet	1	75cm x 90cm	10	2
Hand Towels	4	–		
Leggings (pair)	1	75cm x 114cm		
Perineal Aperture Drape	1	76cm x 170cm		
Instrument Table Cover	1	114cm x 150cm		



9010 Basic Surgical Pack without Mayo


Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Adhesive Towel Drapes	2	100cm x 85cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	



9040 Leggings (pair)


Overall Size	Drapes per Box	Boxes per Case
75cm x 114cm	12	4



1084 Under Buttocks Drape with Pouch

Overall Size	Drapes per Box	Boxes per Case
86cm x 108cm	30	4


- Cuff Size: 25.4cm
- Fluid Collection Pouch



1031 Laparoscopic Drape - Supine Position

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
285cm x 182cm	27cm x 30cm	10	2

- 4 Tube Organisers (ref. 1115)



1032 Laparoscopic Drape - Lithotomy Position


Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
241cm x 228cm	27cm x 30cm	8	2

- 4 Tube Organisers (ref. 1115)
- Perineal Adhesive Aperture: 12cm x 15.5cm

1032P Laparoscopic Drape Pack / Lithotomy Position

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Laparoscopy Drape	1	241cm x 228cm	5	2
Instrument Table Cover		200cm x 150cm		
Hand Towels	4	-		
Mayo Stand Cover	1	80cm x 145cm		

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area




6665 3M™ Ioban™ 2 Abdomino-Perineal Drape - Lithotomy Position

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
200cm x 270cm	35cm x 46cm	10	1

- Tube Organisers
- Perineal Adhesive Aperture: 12cm x.15.5cm

6665P001 Abdomino-Perineal Draping Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Abdomino-Perineal Drape	1	200cm x 270cm	10	1
Instrument Table Cover	1	100cm x 150cm		
Mayo Stand Cover	1	80cm x 145cm		
Towel Drape	1	75cm x 90cm		
Hand Towels	4	-		




2065 3M™ Steri-Drape™ 2 Incise Area Abdomino-Perineal Drape-Lithotomy Position

Overall Size	3M™ Steri-Drape™ 2 Incise Area	Pads per Box	Boxes per Case
195cm x 271cm	41cm x 45cm	8	1

- Tube Organisers
- Perineal Adhesive Aperture: 12cm x 15.5cm

2065P Abdomino-Perineal Draping Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Abdomino-Perineal Drape	1	195cm x 271cm	5	2
Instrument Table Cover	1	100cm x 150cm		
Mayo Stand Cover	1	80cm x 145cm		
Towel Drape	1	75cm x 90cm		
Hand Towels	4	-		



1079 Gynaecological Drapes

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
195cm x 179cm	12cm x 9cm	14	2


- Available as a single box

1079P Gynaecological Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Gynaecology Drapes	1	195cm x 179cm	7	2
Instrument Table Cover	1	100cm x 150cm		
Towel Drapes	1	75cm x 90cm		

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Obstetrics and Gynaecology Drapes




1080 Gynaecological Drape with 1072 Pouch

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
195cm x 179cm	12cm x 9cm	14	2

- Available as a single box

1080P Gynaecological Pack with 1072 Pouch

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Gynaecology Drape	1	195cm x 179cm	7	2
Instrument Table Cover	1	200cm x 150cm		
Towel Drape	1	75cm x 90cm		




6697 Caesarean Section Drape with 3M™ Ioban™ 2 Incise Pouch

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
254cm x 291cm	30cm x 30cm	5	1

- 360-Degree Fluid Collection Pouch
- Tube Organisers

6697P001 Caesarean Section Pack with 3M™ Ioban™ 2 Incise Pouch

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
C-Section Drape	1	254cm x 291cm	10	1
Instrument Table Cover	1	100cm x 150cm		
Mayo Stand Cover	1	80cm x 145cm		
Operation Tape	1	10cm x 55cm		
Hand Towels	4	-		




7966 Caesarean Section Drape with Aperture Pouch

Overall Size	Fenestrated Incise Area	Drapes per Box	Boxes per Case
254cm x 307cm	30cm x 30cm	5	1

- 360-Degree Fluid Collection Pouch
- Tube Organisers

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	




7965 Caesarean Section Drape with 3M™ Steri-Drape™ 2 Incise Pouch

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
254cm x 307cm	30cm x 30cm	5	1

- 360-Degree Fluid Collection Pouch
- Tube Organisers


7965P Caesarean Section Pack with 3M™ Steri-Drape™ 2 Incise Pouch

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
C-Section Drape	1	254cm x 307cm	5	2
Instrument Table Cover	1	100cm x 150cm		
Hand Towels	4	-		
Mayo Stand Cover	1	80cm x 145cm		
Operation Tape	1	10cm x 55cm		



9023 Perineal Aperture Drape


Overall Size	Drapes per Box	Boxes per Case
76cm x 170cm	25	4



6657 Pouch with 3M™ Ioban™ 2 Incise Film

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
87cm x 74cm	30cm x 30cm	10	4


- 360-Degree Fluid Collection Pouch
- Tube Organisers



6658 Large Pouch with 3M™ Ioban™ 2 Incise Film

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
74cm x 74cm	33cm x 43cm	5	4

- Fluid Collection Pouch




6659 Extra Large Pouch with 3M™ Ioban™ 2 Incise Film

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
74cm x 87cm	42cm x 52cm	5	4

- 2 Exit Ports


KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	



1072 Aperture Pouch Drape

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
23cm x 45cm	23cm x 25cm	11cm x 8.5cm	10	4


- Pouch Size: 15.5cm x 23cm



2057 Incise Pouch

Overall Size	3M™ Steri-Drape™ Incise Area	Drapes per Box	Boxes per Case
87cm x 74cm	30cm x 30cm	10	4


- 360-Degree Fluid Collection Pouch
- Tube Organisers



2058 Aperture Pouch

Overall Size	Fenestrated Incise Area	Drapes per Box	Boxes per Case
87cm x 74cm	30cm x 30cm	10	4

- 360-Degree Fluid Collection Pouch
- Tube Organisers



1016 Irrigation Pouch

Overall Size	Drapes per Box	Boxes per Case
50cm x 55cm	10	4

- Adhesive Strip
- 1 Exit Port

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

**Cardiovascular Drapes/
Cardiac Catheterisation Drapes**

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

**Disposable Surgical Gowns
– Standard and Reinforced**

Cross Reference by Product Number

3M™ Urology Drapes

Large amounts of fluids during urological procedures require drapes that can effectively control this fluid and maintain a safe environment for healthcare workers. Drapes that minimise the amount of preparation time before surgery are an advantage as often there is a quick turnaround time between procedures.

Fluid Control Features

- 3M surgical drapes effectively control and contain fluid. This not only reduces potential exposure to infectious body fluids, but keeps the patient dry which can help reduce turnaround time between procedures. Many 3M surgical drapes have fluid collection pouches, and the 3M™ Steri-Drape™ 9000 Range of drapes is absorbent.

3M Range


- 3M surgical drapes are available individually and in specialised packs to cover most surgical procedures. The contents of each draping pack are tailored to match the requirements of each and every procedure. Fully customised packs are also available, to save time and minimise waste.

Impervious to Strikethrough


- A drape that allows fluid transfer, and thus bacteria transfer, compromises the sterile field, which increases the risk of wound infection. That is why a drape that is impervious to liquid strikethrough is so important.

3M™ Ioban™ 2 Antimicrobial Incise Drapes

- Ioban 2 antimicrobial incise drapes contain iodophor which is slowly released onto the skin surface giving continuous and controlled antimicrobial activity throughout the surgical procedure. An Ioban incise area is featured in some speciality drapes. Ioban incise drapes are also available as separate items.




3M™ Urology Drapes




9000 Basic Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Adhesive Towel Drapes	2	100cm x 75cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		




9003 Cystoscopy Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Drape Sheet	1	75cm x 90cm	10	2
Hand Towels	4	–		
Leggings (pair)	1	75cm x 114cm		
Perineal Aperture Drape	1	76cm x 170cm		
Instrument Table Cover	1	114cm x 150cm		



9010 Basic Surgical Pack without Mayo

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Adhesive Towel Drapes	2	100cm x 85cm		
Adhesive Drape Sheet (M)	1	183cm x 183cm		
Adhesive Drape Sheet (L)	1	250cm x 150cm		
Instrument Table Cover	1	200cm x 150cm		




9023 Perineal Aperture Drape

Overall Size	Drapes per Box	Boxes per Case
76cm x 170cm	25	4


KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	

3M™ Urology Drapes



9040 Leggings (pair)

Overall Size	Drapes per Box	Boxes per Case
75cm x 114cm	12	4




1079 Gynaecological Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
195cm x 179cm	12cm x 9cm	14	2

- Available as a single box

1079P Gynaecological Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Gynaecology Drape	1	195cm x 179cm	7	2
Instrument Table Cover	1	100cm x 150cm		
Towel Drape	1	75cm x 90cm		



1080 Gynaecological Drape with 1072 Pouch

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
195cm x 179cm	12cm x 9cm	7	2

- Available as a single box

1080P Gynaecological Pack with 1072 Pouch

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Gynaecology Drape	1	195cm x 179cm	7	2
Instrument Table Cover	1	200cm x 150cm		
Towel Drape	1	75cm x 90cm		



1081 TUR/Urological Drape

Overall Size	Abdominal Adhesive Aperture Size	Drapes per Box	Boxes per Case
152cm x 164cm	6cm	14	2


- Elastic Aperture
- Neoprene Finger Cot
- Large Integrated Fluid Collection Pouch with Filter and Exit Ports

1081P TUR Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
TUR Drape	1	152cm x 164cm	7	2
Instrument Table Cover	1	100cm x 150cm		
Hand Towels	4	-		

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	


3M™ Urology Drapes



1071 Urological Sheet

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
60cm x 90cm	5cm	10	4


- Neoprene Finger Cot with Adhesive
- Adhesive Strip



1072 Aperture Pouch Drape

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
23cm x 45cm	23cm x 25cm	11cm x 8.5cm	10	4

- Pouch Size: 15.5cm x 23cm



2065 Abdomino-Perineal Drape-Lithotomy Position

Overall Size	3M™ Steri-Drape™ 2 Incise Area	Pads per Box	Boxes per Case
195cm x 271cm	41cm x 45cm	8	1

- Tube Organisers
- Perineal Adhesive Aperture: 12cm x 15.5cm

2065P Abdomino-Perineal Draping Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Abdomino-Perineal Drape	1	195cm x 271cm	10	1
Instrument Table Cover	1	100cm x 150cm		
Mayo Stand Cover	1	80cm x 145cm		
Towel Drape	1	75cm x 90cm		
Hand Towels	4	-		



6665 3M™ Ioban™ 2 Abdomino-Perineal Drape - Lithotomy Position

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
200cm x 270cm	35cm x 46cm	10	1


- Tube Organisers
- Perineal Adhesive Aperture: 12cm x 15.5cm

6665P001 Abdomino-Perineal Draping Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Abdomino-Perineal Drape	1	200cm x 270cm	10	1
Instrument Table Cover	1	100cm x 150cm		
Mayo Stand Cover	1	80cm x 145cm		
Towel Drape	1	75cm x 90cm		
Hand Towels	4	-		

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area


3M™ Urology Drapes



1031 Laparoscopic Drape - Supine Position

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
285cm x 182cm	27cm x 30cm	10	2

- Available as a single box
- 4 Tube Organisers (ref. 1115)




1032 Laparoscopic Drape - Lithotomy Position

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
241cm x 228cm	27cm x 30cm	8	2

- 4 Tube Organisers (ref. 1115)
- Perineal Adhesive Aperture: 12cm x 15.5cm
- Available as a single box

1032P Laparoscopic Drape Pack / Lithotomy Position


Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Laparoscopy Drape	1	241cm x 228cm	5	2
Instrument Table Cover		200cm x 150cm		
Hand Towels	4	-		
Mayo Stand Cover	1	80cm x 145cm		



1016 Irrigation Pouch

Overall Size	Units per Box	Boxes per Case
50cm x 55cm	10	4

- Adhesive Strip
- 1 Exit Port



9050 Adhesive Aperture Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
228cm x 300cm	12cm x 16cm	10	2

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

Cardiovascular Drapes/
Cardiac Catheterisation Drapes

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

Disposable Surgical Gowns
– Standard and Reinforced

Cross Reference by Product Number

3M™ Ophthalmic Drapes

Ophthalmic procedures are often carried out under local anaesthesia and therefore the drapes used need to be lightweight and comfortable for the patient. 3M ophthalmic drapes are made from a lightweight and breathable material that allows fluid to run away from the operative site into specifically designed fluid collection pouches, ensuring the drapes remain dry and the patient is comfortable.

Fluid Control Features


- 3M surgical drapes effectively control and contain fluid. This not only reduces potential exposure to infectious body fluids, but keeps the patient dry which can help reduce turnaround time between procedures.

3M Range

- 3M surgical drapes are available individually and in specialised packs to cover most surgical procedures. The contents of each draping pack are tailored to match the requirements of each and every procedure. Fully customised packs are also available, to save time and minimise waste.




3M™ Ophthalmic Drapes



1020 Small Drape with Adhesive Aperture


Overall Size	Adhesive Size	Drapes per Box	Boxes per Case
40cm x 40cm	9.9cm x 12.5cm	10	4

- Circular Aperture Size: 6.3cm



1021 Small Drape with Incise Film


Overall Size	3M™ Steri-Drape™ Incise Area	Drapes per Box	Boxes per Case
40cm x 40cm	9.9cm x 12.5cm	10	4



1022 Refractive Drape

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
40cm x 27cm	7.2cm x 7.2cm	4.4cm x 3.8cm	20	4


- Clear Plastic



1023 Small Drape with Incise Film and Pouch

Overall Size	Steri-Drape incise area	Drapes per Box	Boxes per Case
40cm x 49cm	12.5cm x 10cm	10	4

- 2 Adhesive Patches
- Fluid Collection Pouch




1024 Small Drape with Adhesive Aperture and Pouch

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
39cm x 46cm	12.5cm x 10cm	3.8cm x 6.3cm	10	4

- 2 Adhesive Patches
- Fluid Collection Pouch

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	


3M™ Ophthalmic Drapes



1030 Medium Drape with Adhesive Aperture

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
80cm x 130cm	10cm x 12.5cm	5.1cm x 7.4cm	10	4


- Clear Plastic



1033 Medium Drape with Wide Adhesive Aperture

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
80cm x 130cm	23cm x 11.5cm	15cm x 6cm	10	4


- Clear Plastic



1060 Medium Drape with Incise Film

Overall Size	3M™ Steri-Drape™ Incise Area	Drapes per Box	Boxes per Case
130cm x 130cm	10cm x 12.5cm	10	4

- Clear Plastic




1061 Medium Drape with Incise Film and Pouch

Overall Size	Steri-Drape incise area	Drapes per Box	Boxes per Case
130cm x 130cm	10cm x 12.5cm	10	4

- Fluid Collection Pouch
- Clear Plastic

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	


3M™ Ophthalmic Drapes



1055B Large Ophthalmic Drape with Incise Film and Pouches

Overall Size	3M™ Steri-Drape™ Incise Area	Drapes per Box	Boxes per Case
155cm x 275cm	12cm x 9cm	14	2


- 2 Fluid Collection Pouches
- Face Arch
- Repellent Fabric



1065B Large Ophthalmic Sheet with Aperture and Pouches

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
155cm x 275cm	12cm x 9cm	6cm x 3cm	14	2

- 2 Fluid Collection Pouches
- Face Arch
- Repellent Fabric



1062 Ophthalmic Drape with Incise Film and Pouch

Overall Size	Steri-Drape incise area	Drapes per Box	Boxes per Case
137cm x 137cm	8cm x 10cm	30	2


- Fluid Collection Pouch
- Repellent Fabric

1062P001 Ophthalmology Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Ophthalmology Drape	1	137cm x 137cm	18	1
Instrument Table Cover	1	100cm x 150cm		
Patient Face Protector	1	30cm		

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Ophthalmic Drapes



1064 Ophthalmic Sheet with Aperture and Pouches

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
155cm x 178cm	12cm x 9cm	6cm x 3cm	20	2

- 2 Fluid Collection Pouches
- Repellent Fabric

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

Cardiovascular Drapes/
Cardiac Catheterisation Drapes

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

Disposable Surgical Gowns
– Standard and Reinforced

Cross Reference by Product Number

3M™ Head and Neck Drapes

Head and neck surgery procedures pose a number of challenges for a draping system. Procedures on the ears, nose and throat can be difficult to drape and therefore requires a versatile but comfortable range. Specific needs of some procedures will require speciality drapes, drapes with reliable adhesion so they do not move and, where required, fluid control features.

Superior Strength and Comfort

- The material of the 3M™ Steri-Drape™ 9000 Range of surgical drapes is strong enough to withstand the most arduous procedures, and even the heaviest surgical instruments can be placed on them without concern. They are also soft and pliable for both patient and staff.

3M Range

- 3M surgical drapes are available individually and in specialised packs to cover most surgical procedures. The contents of each draping pack are tailored to match the requirements of each and every procedure. Fully customised packs are also available, to save time and minimise waste.

Reliable Adhesion



- 3M is a leader in adhesive technology, and has applied this expertise to provide surgical drapes with consistently reliable adhesion to the skin or wound edge. Even during heavy retraction and irrigation, 3M incise drapes adhere securely to wound edges, providing a barrier to bacteria on the skin.

Fluid Control Features


- 3M surgical drapes effectively control and contain fluid. This not only reduces potential exposure to infectious body fluids, but keeps the patient dry which can help reduce turnaround time between procedures. Many 3M surgical drapes have fluid collection pouches, and the Steri-Drape 9000 range of drapes is absorbent.

3M™ Ioban™ 2 Antimicrobial Incise Drapes

- Ioban 2 antimicrobial incise drapes contain iodophor which is slowly released onto the skin surface giving continuous and controlled antimicrobial activity throughout the surgical procedure. An Ioban incise area is featured in some speciality drapes. Ioban incise drapes are also available as separate items.




3M™ Head and Neck Drapes




9007 ENT Surgical Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tapes	2	10cm x 55cm	8	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Drape Sheet	1	75cm x 90cm		
Adhesive Split Sheet	1	150cm x 250cm		
Adhesive Towel Drape	1	100cm x 75cm		
Instrument Table Cover	1	150cm x 150cm		



9008 Head and Neck Pack

Product	Quantity	Overall Size	Packs per Box	Boxes per Case
Operation Tape	1	10cm x 55cm	7	2
Hand Towels	4	–		
Mayo Stand Cover	1	80cm x 145cm		
Headwrap	1	114cm x 90cm		
Adhesive Split Sheet	1	150cm x 250cm		
Instrument Table Cover	1	200cm x 150cm		




9051 Thyroid Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
183cm x 315cm	12cm x 12cm	15	2


KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	

3M™ Head and Neck Drapes



6687 Craniotomy Drape with 3M™ Ioban™ 2 Incise Film


Overall Size	Incise Area	Drapes per Box	Boxes per Case
231cm x 432cm	36.8cm x 21.4cm	10	1



9047 Neuro Split Sheet with Pouch


Overall Size	Adhesive Split Size	Drapes per Box	Boxes per Case
183cm x 152cm	16cm x 40cm	15	2

- Fluid Collection Pouch
- 2 Exit Ports



9048 Split Sheet

Overall Size	Adhesive Split Size	Drapes per Box	Boxes per Case
228cm x 260cm	20cm x 40cm	15	2



The 9047 and the 9048 (Split Sheet) will together create a Neuro Drape

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

**Cardiovascular Drapes/
Cardiac Catheterisation Drapes**

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

**Disposable Surgical Gowns
– Standard and Reinforced**

Cross Reference by Product Number

3M™ Neurosurgery Drapes

Neurosurgical procedures are often long and involve large amounts of blood and irrigation fluids. Drapes that can adequately collect this fluid and maintain a safe environment for healthcare workers are an advantage and will ensure the risk of infection is also reduced. Ease of application is also imperative as these procedures are often emergencies.

Fluid Control Features

- 3M surgical drapes effectively control and contain fluid. This not only reduces potential exposure to infectious body fluids, but keeps the patient dry which can help reduce turnaround time between procedures.

Superior Strength and Comfort

- The material of the 3M™ Steri-Drape™ 9000 Range of surgical drapes is strong enough to withstand the most arduous procedures, and even the heaviest surgical instruments can be placed on them without concern. They are also soft and pliable for both patient and staff.

Impervious to Strikethrough


- A drape that allows fluid transfer, and thus bacteria transfer, compromises the sterile field, which increases the risk of wound infection. That is why a drape that is impervious to liquid strikethrough is so important.

3M Range


- 3M surgical drapes are available individually and in specialised packs to cover most surgical procedures. The contents of each draping pack are tailored to match the requirements of each and every procedure. Fully customised packs are also available, to save time and minimise waste.

3M™ Ioban™ 2 Antimicrobial Incise Drapes

- Ioban 2 antimicrobial incise drapes contain iodophor which is slowly released onto the skin surface giving continuous and controlled antimicrobial activity throughout the surgical procedure. An Ioban incise area is featured in some speciality drapes. Ioban incise drapes are also available as separate items.




3M™ Neurosurgery Drapes



6687 Craniotomy Drape with 3M™ Ioban™ 2 Incise Film


Overall Size	Incise Area	Drapes per Box	Boxes per Case
231cm x 432cm	36.8cm x 21.4cm	10	1



9047 Neuro Split Sheet with Pouch


Overall Size	Adhesive Split Size	Drapes per Box	Boxes per Case
183cm x 152cm	16cm x 40cm	15	2

- Fluid Collection Pouch
- 2 Exit Ports



9048 Split Sheet

Overall Size	Adhesive Split Size	Drapes per Box	Boxes per Case
228cm x 260cm	20cm x 40cm	15	2



The 9047 and the 9048 (Split Sheet) will together create a Neuro Drape

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

**Cardiovascular Drapes/
Cardiac Catheterisation Drapes**

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

**Disposable Surgical Gowns
– Standard and Reinforced**

Cross Reference by Product Number


Accessories and other Drapes

Whatever your draping requirements, 3M can supply an additional range of accessories to complement the speciality drapes and drape packs available.

The range of adhesive and non-adhesive individually packed drapes combine the benefits of absorbent, impervious fabric and the flexibility of cost effective draping options. In addition, the range of plastic drapes and equipment covers ensure all surgical procedural needs are met.



3M™ Fabric Drapes



Adhesive Split Sheets

Split Size 7 x 70cm


9043

9044

Split Size 20 x 100cm

9045

Overall Size	Drapes per Box	Boxes per Case
150cm x 183cm	30	2
150cm x 250cm	25	2
228cm x 260cm	15	2



9053 Extremity Drape

Overall Size

228cm x 300cm

Aperture Size


6cm x 8cm

Drapes per Box

10

Boxes per Case

2



9050 Adhesive Aperture Drape

Overall Size

228cm x 300cm

Drapes per Box


10

Boxes per Case

2


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Fabric Drapes




9046 Bilateral Split Sheet

Overall Size	Aperture Split	Drapes per Box	Boxes per Case
183cm x 200cm	15cm x 35cm	15	2




9051 Thyroid Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
183cm x 315cm	12cm x 12cm	15	2



9026 Aperture Drape

Overall Size	Drapes per Box	Boxes per Case
55cm x 55cm	100	4




9029 Adhesive Aperture Drape

Overall Size	Drapes per Box	Boxes per Case
75cm x 75cm	40	4


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Fabric Drapes




Adhesive Aperture Drapes

	Overall Size	Drapes per Box	Boxes per Case
Centred Aperture 9028	75cm x 75cm	40	4
9030	120cm x 150cm	20	4
Off-Centre Aperture 9031	150cm x 183cm	25	2




Adhesive Towel Drapes

	Overall Size	Drapes per Box	Boxes per Case
9084	75cm x 75cm	40	4
9085	100cm x 85cm	30	4
9086	100cm x 75cm	40	4



Adhesive Drape Sheets

	Overall Size	Drapes per Box	Boxes per Case
9087	183cm x 183cm	25	2
9088	250cm x 150cm	25	2
9090	300cm x 228cm	12	2




9052 Hand Drape

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
360cm x 200cm	3cm	9	2


KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	

3M™ Fabric Drapes




9054 Foot / Hand Drape

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
200cm x 300cm	3cm	12	2




Non-Adhesive Drape Sheets/ Instrument Table Covers

	Overall Size	Drapes per Box	Boxes per Case
9067	37cm x 45cm	150	4
9069	75cm x 90cm	30	4
9072	100cm x 150cm	30	2
9075	150cm x 170cm	15	2
9077	200cm x 150cm	15	2
9079	200cm x 228cm	20	2



9035 Mayo Stand Cover

Overall Size	Drapes per Box	Boxes per Case
80cm x 145cm	40	2



9040 Leggings (pair)

Overall Size	Drapes per Box	Boxes per Case
75cm x 114cm	12	4


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Fabric Drapes



Leggings


	Overall Size	Drapes per Box	Boxes per Case
9041	35cm x 60cm	50	4
9042	35cm x 100cm	35	4



9062 Laparotomy Drape (Transverse)


Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
254cm x 285cm	30cm x 10cm	12	2

- Absorbent Impervious Material




9063 Laparotomy Drape

Overall Size	Adhesive Aperture Size	Drapes per Box	Boxes per Case
254cm x 285cm	10cm x 30cm	10	2



9099 Operation Tape

Overall Size	Items per Box	Boxes per Case
10cm x 55cm	250	2




9095 Absorbent Pad

Overall Size	Items per Box	Boxes per Case
70cm x 50cm	50	4


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Clear Plastic Drapes




1000 Small Towel Drape

Overall Size	Adhesive Strip Width	Drapes per Box	Boxes per Case
45cm x 30cm	3cm	10	4



1010 Large Towel Drape


Overall Size	Adhesive Strip Width	Drapes per Box	Boxes per Case
45cm x 60cm	3cm	10	4



1020 Small Drape with Adhesive Aperture


Overall Size	Adhesive Size	Drapes per Box	Boxes per Case
40cm x 40cm	9.9cm x 12.5cm	10	4

• Aperture Size: 6.3cm




1021 Small Drape with Incise Film

Overall Size	3M™ Steri-Drape™ Incise Area	Drapes per Box	Boxes per Case
40cm x 40cm	9.9cm x 12.5cm	10	4



1030 Medium Drape with Adhesive Aperture

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
80cm x 130cm	10cm x 12.5cm	5.1cm x 7.4cm	10	4




1033 Medium Drape with Wide Adhesive Aperture

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
80cm x 130cm	23cm x 11.5cm	15cm x 6cm	10	4


- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area

3M™ Clear Plastic Drapes



1060 Medium Drape with Incise Film


Overall Size	3M™ Steri-Drape™ Incise Area	Drapes per Box	Boxes per Case
130cm x 130cm	10cm x 12.5cm	10	4



1071 Urological Sheet

Overall Size	Aperture Size	Drapes per Box	Boxes per Case
60cm x 90cm	5cm	10	4


- Neoprene Finger Cot with Adhesive
- Adhesive Strip



1015 U-Drape

Overall Size	Drapes per Box	Boxes per Case
120cm x 130cm	10	4


- U-Slot Aperture with Adhesive



1019 Long U-Drape

Overall Size	Drapes per Box	Boxes per Case
120cm x 180cm	10	4

- U-Slot Aperture with Adhesive
- Adhesive Patch



Wound Edge Protector

	Overall Size	Ring Diameter	Packs per Box	Boxes per Case
1073	90cm x 90cm	12.5cm	10	4
1074	90cm x 90cm	17.5cm	10	4
1075	90cm x 90cm	22.5cm	10	4
1076	90cm x 90cm	27.5cm	10	4
1077	90cm x 90cm	9.5cm	10	4

- 4 Adhesive Patches
- 1 Flexible Ring

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	


Accessories



1003 Isolation Bag


Overall Size	Drapes per Box	Boxes per Case
49cm x 49cm	10	4

- Drawstring Closure
- Clear Plastic




1115 Tube Organiser

Overall Size	Units per Box	Boxes per Case
14cm x 4cm	100	4



9097 Instrument Pouch (L)


Overall Size	Drapes per Box	Boxes per Case
33cm x 38cm	100	2



1018 Instrument Pouch

Overall Size	Drapes per Box	Boxes per Case
17cm x 29cm	10	4


- 2 Adhesive Strips
- 2 Compartments



1018L Instrument Pouch


Overall Size	Drapes per Box	Boxes per Case
24cm x 46cm	10	4

- 2 Adhesive Strips
- 3 Compartments



9096 Arm Cover

Overall Size	Items per Box	Boxes per Case
54cm	100	1




Stockinette

	Overall Size	Drapes per Box	Boxes per Case
9109	23cm x 71cm	20	1
9110	30cm x 122cm	20	1

- KEY**
- Absorbent Impervious Fabric
 - Repellent Fabric
 - Clear Plastic Film
 - Blue Plastic Film
 - Extra Absorbent Area
 - 3M™ Ioban™ 2 Antimicrobial Incise Film
 - 3M™ Steri-Drape™ Incise Film
 - Clear Elastomeric Film
 - Adhesive Area


3M™ Fluid Collection Pouches



1007 Arthroscopy Fluid Collection Pouch

Overall Size	Drapes per Box	Boxes per Case
115cm x 58cm	10	4


- Adhesive Strips
- 1 Exit Port



1084 Under Buttocks Drape with Pouch

Overall Size	Drapes per Box	Boxes per Case
86cm x 108cm	30	4


- Fluid Collection Pouch



1027 Irrigation Pouch

Overall Size	Drapes per Box	Boxes per Case
45cm x 30cm	10	4


- Adhesive Strip
- Exists Non-Sterile (ref. 1026)



2057 Incise Pouch

Overall Size	3M™ Steri-Drape™ Incise Area	Drapes per Box	Boxes per Case
87cm x 74cm	30cm x 30cm	10	4


- 360-Degree Fluid Collection Pouch
- Tube Organisers



6657 Pouch with 3M™ Ioban™ 2 Incise Film

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
87cm x 74cm	30cm x 30cm	10	4


- 360-Degree Fluid Collection Pouch
- Tube Organisers



6658 Large Pouch with Ioban 2 incise film

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
74cm x 74cm	33cm x 43cm	5	4

- Fluid Collection Pouch




6659 Extra Large Pouch with Ioban 2 incise film

Overall Size	Ioban 2 incise area	Drapes per Box	Boxes per Case
74cm x 87cm	42cm x 52cm	5	4

- 2 Exit Ports

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	


3M™ Fluid Collection Pouches



1067 U-Pouch

Overall Size	Drapes per Box	Boxes per Case
87cm x 74cm	5	4


- U-Slot Aperture with Adhesive
- 2 Exit Ports



1016 Irrigation Pouch

Overall Size	Drapes per Box	Boxes per Case
50cm x 55cm	10	4


- Adhesive Strip
- 1 Exit Port



1072 Aperture Pouch Drape

Overall Size	Adhesive Size	Oval Aperture Size	Drapes per Box	Boxes per Case
23cm x 45cm	23cm x 25cm	11cm x 8.5cm	10	4

- Pouch Size: 15.5cm x 23cm



2058 Aperture Pouch

Overall Size	Fenestrated Incise Area	Drapes per Box	Boxes per Case
87cm x 74cm	30cm x 30cm	10	4

- 360-Degree Fluid Collection Pouch
- Tube Organisers


3M™ Equipment Cover Drapes



1100E Camera/Laser Drape

Overall Size	Drapes per Box	Boxes per Case
250cm x 13cm	50	1


- With Elastoflex™ KG Semi-Opaque Fenestration



1012 Fluoroscope Drape

Overall Size	Drapes per Box	Boxes per Case
90cm x 110cm	10	4


- 2 Adhesive Strips
- 2 Adhesive Patches



1101 One Piece C-Arm Drape

Overall Size	Drapes per Box	Boxes per Case
265cm x 150cm	10	2


- Hook and Loop Fasteners, Poly Straps and Adhesive Strips



X-ray Cassette Drape

	Overall Size	Drapes per Box	Boxes per Case
1008	36cm x 49cm	25	1
1009	56cm x 59cm	25	1
1011	57.5cm x 62.5cm	25	2

- Adhesive Closure




1013 X-ray Image Intensifier Drape

Overall Size	Drapes per Box	Boxes per Case
165cm x 106cm	10	4

- 4 Adhesive Patches

KEY	Absorbent Impervious Fabric	Repellent Fabric	Clear Plastic Film	Blue Plastic Film	Extra Absorbent Area
	3M™ Ioban™ 2 Antimicrobial Incise Film	3M™ Steri-Drape™ Incise Film	Clear Elastomeric Film	Adhesive Area	

3M™ Equipment Cover Drapes



Banded Bag

	Overall Size	Drapes per Box	Boxes per Case
1096	50cm	25	1
1097	66cm	25	1
1098	100cm	25	1
1099	127cm	25	1

• Latex-free Elasticated Opening

KEY

-  Absorbent Impervious Fabric
-  Repellent Fabric
-  Clear Plastic Film
-  Blue Plastic Film
-  Extra Absorbent Area
-  3M™ Ioban™ 2 Antimicrobial Incise Film
-  3M™ Steri-Drape™ Incise Film
-  Clear Elastomeric Film
-  Adhesive Area

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

Cardiovascular Drapes/
Cardiac Catheterisation Drapes

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

**Disposable Surgical Gowns
– Standard and Reinforced**

Cross Reference by Product Number

3M™ High Performance Disposable Surgical Gowns – Standard and Reinforced

Non-woven and made of four layers for optimal patient and wearer protection.

- Lightweight and breathable for maximum comfort.
- High barrier performance.
- Good strength characteristics.
- Good breathability.
- Low linting.




3M™ Disposable Surgical Gowns – Standard and Reinforced

Traditional 3M disposable surgical gown.

- Light and soft material keeps the wearer dry and cool.
- Standard and reinforced designs.




3M™ High Performance Disposable Surgical Gowns
 – Standard and Reinforced



**High Performance Disposable Surgical Gown
 – Standard**

	Size	Boxes per Case	Items per Box
7691C	Medium	1	28
7692C	Large	1	26
7693C	XLarge	1	24
7694C	XXLarge	1	22




**High Performance Disposable Surgical Gown
 – Reinforced**

	Size	Boxes per Case	Items per Box
7695C	Medium	1	22
7696C	Large	1	20
7697C	XLarge	1	18
7698C	XXLarge	1	16


KEY  Impervious Fabric  Repellent Fabric

3M™ Disposable Surgical Gowns – Standard and Reinforced



Disposable Surgical Gown

	Size	Boxes per Case	Items per Box
7591	Medium	1	25
7592	Large	1	25
7593	XLarge	1	25
7594	XXLarge	1	25



Reinforced Disposable Surgical Gown

	Size	Boxes per Case	Items per Box
7595	Medium	1	25
7596	Large	1	23
7597	XLarge	1	23
7598	XXLarge	1	23

KEY  Impervious Fabric  Repellent Fabric

Custom Procedure Trays

Incise Drapes

General Surgery Drapes

Orthopaedic Drapes

**Cardiovascular Drapes/
Cardiac Catheterisation Drapes**

Obstetrics and Gynaecology Drapes

Urology Drapes

Ophthalmic Drapes

Head and Neck Drapes

Neurosurgery Drapes

Accessories and other Drapes

**Disposable Surgical Gowns
– Standard and Reinforced**

Cross Reference by Product Number

3M™ Custom Procedure Trays		Page
-	Custom Procedure Trays	10

3M™ Incise Drapes		Page	3M™ Incise Drapes		Page
1035	Incise Drape	15	2051	Incise Drape	16
1037	Incise Drape	15	6635	Antimicrobial Incise Drape	13
1040	Incise Drape	15	6640	Antimicrobial Incise Drape	13
1050	Incise Drape	15	6648	Antimicrobial Incise Drape	13
1051	Incise Drape	15	6650	Antimicrobial Incise Drape	13
2035	Incise Drape	16	6651	Antimicrobial Incise Drape	13
2037	Incise Drape	16	6640EZ	Antimicrobial Incise Drape	14
2040	Incise Drape	16	6648EZ	Antimicrobial Incise Drape	14
2045	Incise Drape	16	6650EZ	Antimicrobial Incise Drape	14
2050	Incise Drape	16	6651EZ	Antimicrobial Incise Drape	14

3M™ Speciality Packs/Drapes		Page	3M™ Speciality Packs/Drapes		Page
1000	Small Towel Drape	76	1062P001	Ophthalmology Pack	59
1003	Isolation Bag	24, 78	1064	Ophthalmic Sheet with Aperture and Pouches	60
1007	Arthroscopy Fluid Collection Pouch	36, 79	1065B	Large Ophthalmic Sheet with Aperture and Pouches	59
1008	X-ray Cassette Drape	81	1067	U-Pouch	36, 80
1009	X-ray Cassette Drape	81	1071	Urological Sheet	53, 77
1010	Large Towel Drape	76	1072	Aperture Pouch Drape	48, 53, 80
1011	X-ray Cassette Drape	81	1073	Wound Edge Protector	24, 77
1012	Fluoroscope Drape	81	1074	Wound Edge Protector	24, 77
1013	X-ray Image Intensifier Drape	81	1075	Wound Edge Protector	24, 77
1014	Patient Isolation Drape	33	1076	Wound Edge Protector	24, 77
1015	U-Drape	34, 77	1077	Wound Edge Protector	24, 77
1016	Irrigation Pouch	23, 36, 48, 54, 80	1079	Gynaecological Drape	45, 52
1017	Isolation Drape with Incise Film and Pouch	33	1079P	Gynaecological Pack	45, 52
1017P001	Patient Isolation Pack	33	1080	Gynaecological Drape with 1072 Pouch	46, 52
1018	Instrument Pouch	78	1080P	Gynaecological Pack with 1072 Pouch	46, 52
1018L	Instrument Pouch – Long	78	1081	TUR/Urological Drape	52
1019	Long U-Drape	34, 77	1081P	TUR Pack	52
1020	Small Drape with Adhesive Aperture	57, 76	1084	Under Buttocks Drape with Pouch	23, 44, 79
1021	Small Drape with Incise Film	57, 76	1090	Femoral Angiography Drape	39
1022	Refractive Drape	57	1092	Minor Procedure Drape	22
1023	Small Drape with Incise Film and Pouch	57	1094	Arthroscopy Sheet with Pouch	35
1024	Small Drape with Adhesive Aperture and Pouch	57	1096	Banded Bag	82
1027	Irrigation Pouch	79	1097	Banded Bag	82
1030	Medium Drape with Adhesive Aperture	58, 76	1098	Banded Bag	82
1031	Laparoscopic Drape – Supine Position	26, 44, 54	1099	Banded Bag	82
1032	Laparoscopic Drape – Lithotomy Position	26, 44, 54	1100E	Camera/Laser Drape	81
1032P	Laparoscopic Drape Pack / Lithotomy Position	26, 44, 54	1101	One Piece C-Arm Drape	81
1033	Medium Drape with Wide Adhesive Aperture	58, 76	1115	Tube Organiser	23, 36, 78
1055B	Large Ophthalmic Drape with Incise Film and Pouches	59	1191	Femoral Angiography Drape	40
1060	Medium Drape with Incise Film	58, 77	1194	Arthroscopy Sheet with Pouch	35
1061	Medium Drape with Incise Film and Pouch	58	1195	Arthroscopy Pack	35
1062	Ophthalmic Drape with Incise Film and Pouch	59	2057	Incise Pouch	23, 48, 79

3M™ Speciality Packs/Drapes		Page	3M™ Speciality Packs/Drapes		Page
2058	Aperture Pouch	23, 48, 80	9029	Adhesive Aperture Drape	72
2065	Abdomino-Perineal Drape – Lithotomy Position	25, 45, 53	9030	Adhesive Aperture Drape	73
2065P	Abdomino-Perineal Draping Pack	25, 45, 53	9031	Adhesive Aperture Drape	73
6617	Isolation Drape with 3M™ Ioban™ 2 Film and Pouch	34	9035	Mayo Stand Cover	74
6617P001	Patient Isolation Pack with 3M™ Ioban™ 2 Film and Pouch	34	9040	Leggings (Pair)	21, 44, 52, 74
6619	Large Isolation Drape with 3M™ Ioban™ 2 Film and Pouch	34	9041	Leggings	33, 75
6657	Pouch with 3M™ Ioban™ 2 Incise Film	24, 47, 79	9042	Leggings	33, 75
6658	Large Pouch with 3M™ Ioban™ 2 Incise Film	24, 47, 79	9043	Adhesive Split Sheet	31, 71
6659	Extra Large Pouch with 3M™ Ioban™ 2 Incise Film	24, 47, 79	9044	Adhesive Split Sheet	31, 71
6665	3M™ Ioban™ 2 Abdomino-Perineal Drape – Lithotomy Position	26, 45, 53	9045	Adhesive Split Sheet	31, 71
6665P001	Abdomino-Perineal Draping Pack	26, 45, 53	9046	Bilateral Split Sheet	39, 72
6677	Cardiovascular Sheet with 3M™ Ioban™ 2 Incise Film	39	9047	Neuro Split Sheet with Pouch	64, 67
6681	3M™ Ioban™ 2 Cardiovascular Drape	40	9048	Split Sheet	64, 67
6687	Craniotomy Drape with 3M™ Ioban™ 2 Incise Film	64, 67	9050	Adhesive Aperture Drape	21, 54, 71
6697	Caesarean Section Sheet with 3M™ Ioban™ 2 Incise Pouch	46	9051	Thyroid Drape	21, 63, 72
6697P001	Caesarean Section Pack with 3M™ Ioban™ 2 Incise Pouch	46	9052	Hand Drape	25, 31, 73
7965	Caesarean Section Drape with 3M™ Steri-Drape™ 2 Incise Pouch	47	9053	Extremity Drape	32, 71
7965P	Caesarean Section Pack with 3M™ Steri-Drape™ 2 Incise Pouch	47	9054	Foot/Hand Drape	25, 32, 74
7966	Caesarean Section Sheet with Aperture Pouch	46	9055	One-Piece Hip Drape with Side Pouches	32
9000	Basic Surgical Pack	19, 29, 43, 51	9062	Laparotomy Drape (Transverse)	22, 75
9000L	Basic Surgical Pack	19	9063	Laparotomy Drape	22, 33, 75
9001	Paediatric Pack	21	9067	Drape Sheet/Instrument Table Cover	74
9002	Orthopaedic Surgical Pack without Mayo	29	9069	Drape Sheet/Instrument Table Cover	74
9003	Cystoscopy Surgical Pack	43, 51	9072	Drape Sheet/Instrument Table Cover	74
9004	Extremity Surgical Pack	29	9075	Drape Sheet/Instrument Table Cover	74
9007	ENT Surgical Pack	19, 39, 63	9077	Drape Sheet/Instrument Table Cover	74
9008	Head and Neck Pack	20, 63	9079	Drape Sheet/Instrument Table Cover	74
9010	Basic Surgical Pack without Mayo	20, 31, 43, 51	9084	Adhesive Towel Drape	73
9012	Orthopaedic Surgical Pack	30	9085	Adhesive Towel Drape	73
9012L	Orthopaedic Surgical Pack	30	9086	Adhesive Towel Drape	73
9014	Extremity Surgical Pack	30	9087	Adhesive Drape Sheet	73
9015	Hand Pack	22, 31	9088	Adhesive Drape Sheet	73
9018	Surgical Split Sheet Pack	20, 40	9090	Adhesive Drape Sheet	73
9023	Perineal Aperture Drape	47, 51	9095	Absorbent Pad	75
9026	Aperture Drape	72	9096	Arm Cover	78
9028	Adhesive Aperture Drape	73	9097	Instrument Pouch (L)	78
			9099	Operation Tape	75
			9109	Stockinette	36, 78
			9110	Stockinette	36, 78
			9196	Shoulder Split Sheet with Pouch	35

Disposable Surgical Gowns		Page	Disposable Surgical Gowns		Page
7591	Standard Disposable Surgical Gown	86	7691C	High Performance Disposable Surgical Gown – Standard	85
7592	Standard Disposable Surgical Gown	86	7692C	High Performance Disposable Surgical Gown – Standard	85
7593	Standard Disposable Surgical Gown	86	7693C	High Performance Disposable Surgical Gown – Standard	85
7594	Standard Disposable Surgical Gown	86	7694C	High Performance Disposable Surgical Gown – Standard	85
7595	Reinforced Disposable Surgical Gown	86	7695C	High Performance Disposable Surgical Gown – Reinforced	85
7596	Reinforced Disposable Surgical Gown	86	7696C	High Performance Disposable Surgical Gown – Reinforced	85
7597	Reinforced Disposable Surgical Gown	86	7697C	High Performance Disposable Surgical Gown – Reinforced	85
7598	Reinforced Disposable Surgical Gown	86	7698C	High Performance Disposable Surgical Gown – Reinforced	85

3M Health Care

Unique Drape Exchange Scheme

3M Health Care recognises that lists and procedures are subject to last minute changes. Therefore your 3M representative will be happy to exchange any uncontaminated opened and unused 3M surgical drape or surgical gown, free of charge, for a new sterile drape or gown.

For further information, samples, or to arrange a free demonstration of gowns and draping systems, contact Customer Services on 01509 613151, or email infectionprevention.uk@mmm.com.

3M Health Care not only provides quality products but also appropriate support services which enable correct product use and optimised outcomes.

3M strives to contribute to one of the ultimate goals of all health care organisations: infection free patients and health care professionals, throughout the continuum of care.

For further information on 3M products please visit our website:
www.3mhealthcare.co.uk



Asepsis Products
3M Health Care Limited

3M House
Morley Street
Loughborough
Leicestershire LE11 1EP
Tel: (01509) 611611
Fax: (01509) 237288

3M Ireland

3M House, Adelphi Centre
Dun Laoghaire
Co. Dublin, Ireland
Tel: (01) 280 3555
Fax: (01) 280 3509

3M, Coban, Steri-Drape and Ioban
are trademarks of the 3M Company.
©3M Health Care Limited 2006.

E-mail: infectionprevention.uk@mmm.com