

Colour versions of upholstery



P – orange

M – corn yellow



B – brilliant blue

F – ocean green



Height adjustment of the chair

Tilting of the back rest

Height phases

Tilting phases

Technical parameters

Dimensions	
Total length:	1920 mm
Total width:	860 mm
Electric adjustment	
Back rest adjustment angle:	-12° / +70°
Seat adjustment:	+5° / +8° (automatic)
Minimum height:	630 mm
Maximum height:	1175 mm
Maximum load	
Maximum patient weight:	175 kg

Foot rest increases comfort of mother



VITA

BIRTHING AND EXAMINATION CHAIR



BORCAD cz s.r.o.

CZ – 739 45 Fryčovice 673
Tel.: +420 338 640 831
Fax: +420 558 658 087

email: borcad@borcad.cz
www.borcad.cz

BORCAD®

healthy
medical
www.hmed.com.hk
Sales Hotline: 8338 3333

VITA

BIRTHING AND EXAMINATION CHAIR



Comfortable access to the mother
Friendly design with compact shapes
Quick and intuitive control
Very easy cleaning

BORCAD®

Advantages

Very easy maintenance and cleaning




Easy cleaning

Reduction of time necessary for cleaning thanks to seamless upholstery and compact, easily accessible areas.



Leg holder

Adjusting the leg holder horizontally as well as vertically using one controller and one hand.



Accessory with smart storage

Birthing and examination positions

Easy access to the mother from all the sides is one of the main benefits of the VITA chair.



A very low access-height of the chair enables easy access to the birthing chair also for a small woman in childbirth.



Accessories



PKAZP001 X Foot rest - removable



PKAZP005 Baby plate - removable



ZP-Z5 X Fixation strap (pair)



PKAZP008 X Side supports (pair)



SIAZP003 Telescopic infusion stand



PKAZP012 X wedged headrest



PKAZP017 X Head rest - semicircular



PKAZP003 X Head rest



PKAZP004 Holder of infusion stand



ZK05 X Doctor's mobile chair

PKAZP016 equipotential coupling