




Baby scale "Easy" with fold function

8320

- ▶ Foldable baby scale with integrated carrying handles
- ▶ Lightweight and handy
- ▶ Eight-point measuring system guarantees exact weight evaluation
- ▶ Eight load cells
- ▶ Easy and fast opening: Pull the handles apart and remove the base



Features

The first baby scale with a handy fold function for parents, midwives and pharmacies can be reduced easily to the half of its length and is impressive because of its low weight. Furthermore the scale is extra handy and easy to transport thanks to the integrated carry handle and a beauty case-shaped design. Even when the child is very fitful the eight-point-weighing guarantees an every time accurate determination of weight which is completed by an easy to read display. The baby scale has an environment-friendly auto-off function for a long battery life.

- ▶ Approved to class III, Medical device directive (MPG) class 1
- ▶ Attractive baby scale with very fine division 2 g / 5 g
- ▶ Large LCD display, 20 mm digit height
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function: The weight value will be held for a few seconds after unloading
- ▶ Baby tray and housing of one piece
- ▶ Analogue and digital measurements with baby height rod after possible
- ▶ Units in kg or lb
- ▶ Energy-saving auto-off function
- ▶ Mains or battery operation

Products

Ordernumber Variation	Approvable	Division g	Maximum load kg
8320.01.001 Baby scale "Easy" with fold function	No	10	15
