



# Baby scale "Home"

8310.01

- ▶ Maximum capacity: 20 kg
- ▶ Division: 10 g
- ▶ Hold function
- ▶ Energy-saving auto-off function

## Features

---

- ▶ Attractive baby scale with very fine division 2 g / 5 g
- ▶ Large LCD display, 20 mm digit height
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function: The weight value will be held for a few seconds after unloading
- ▶ Baby tray and housing of one piece
- ▶ Analogue and digital measurements with baby height rod after possible
- ▶ Units in kg or lb
- ▶ Energy-saving auto-off function
- ▶ Mains or battery operation

## Products

---

| Ordernumber<br>Variation         | Approvable | Division<br>g | Maximum load<br>kg |
|----------------------------------|------------|---------------|--------------------|
| 8310.01.001<br>Baby scale "Home" | No         | 10            | 20                 |

---