



# Baby scale "Multina Comfort"

8352

- ▶ Detachable baby tray
- ▶ With attachable baby height rod (optional)
- ▶ Simple operation with only one button
- ▶ Scale can be used to weigh children up to 20 kg



## Features

---

- ▶ Simple operation with only one button
- ▶ With detachable baby tray
- ▶ Scale can be used for small children up to 20 kg
- ▶ Easy to read LCD-display, 13 mm digit height
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function: The weight value will be held for a few seconds after unloading
- ▶ Units in kg or lb
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 9 V battery (included in scope of delivery)

## Products

---

Ordernumber Variation	Approvable	Division g	Maximum load kg
8352.01.001 Baby scale "Multina Comfort"	No	10	20

---