



# Fitness scale

7850

- ▶ Measurement of muscle, fat and water components, BMI and caloric needs per day
- ▶ Printout of measurement values
- ▶ Simple self-explanatory menu guide on the glass surface
- ▶ Clothing weight value seasonally adjustable



SOEHNLE PROFESSIONAL	
24-03-2010 09:50	
Height:	160cm
Age:	26
Sex:	female
Weight with clothing:	55.4kg
Weight without clothing:	52.9kg
BMI:Gewicht ok	20.7

## Features

---

- ▶ Measurement of muscle, fat and water components in the body
- ▶ Impedance measurement via the hands, no need to take off the shoes
- ▶ Optimized electronics with additional settings
- ▶ Large LCD display
- ▶ Indication of time in standby mode
- ▶ Simple self-explanatory menu guide on the glass surface
- ▶ No limit on measurement for over- and underweight people
- ▶ Clothing weight value seasonally adjustable
- ▶ Automatic data printout of the analysis measurements
- ▶ Company information on printout individually adjustable with included software in the package
- ▶ Extra large weighing platform with improved surface
- ▶ No assembly required - the scale is immediately ready for use
- ▶ Energy saving auto-off function
- ▶ Switchable to kg / st / lb units
- ▶ Mains operation

Available languages:

German, English, French, Spanish, Italian, Greek, Turkish, Hungarian, Dutch, Czech, Danish, Finnish, Swedish, Norwegian, Slovak, Polish, Russian

## Products

---

Ordernumber Variation	Approvable	Division g	Maximum load kg	Platform size mm
7850.30.001 Fitness scale	No	100	180	365 x 490

---