



Physician Scale

H151-00-5 (Class III); H150-00-5 (Non-Class III)

Lightweight and slender, the **Healthweigh®** Physician Scale is a professional tool providing precision weighing information. The design patented scale is multifunctional: height and weight are indicated in a single task operation for BMI calculation. The adjustable height rod is suitable for children and adults.

Features

- Movement compensation technology
- BMI function
- Zero/Hold/Print/kg/lb
- 5 digit 3/4" LCD/Positive FSTN
- Automatic Reweigh
- Updatable firmware
- Height Range: 70 cm – 204 cm
- Ports: RS232 and USB
- Optional: Bluetooth

Power

- 230v/120v input, DC 9V output
- Auto power off
- 6 X 1.5AA batteries
- Low battery indicator

Technical Data

Class III & Non-Class III	
Capacity	250kg / 550lb
Graduation	0.1kg / 0.2lb
Total Weight	12kg / 26.5lb
Product Size	567 x 375 x 1331 mm 22.32" x 14.76" x 52.40"
Platform Size	375 x 375 x 78 mm 14.76" x 14.76" x 3"

Class III EMR
READY

An outstanding platform

Aluminum cast coat makes the scale safe and robust. The platform is wide and coated with silicone mat for patient comfort.



reddot design award
winner 2011



SHEKEL All the weigh!