



Height rod, foldable

5003

- ▶ Handy height rod for mobile use
- ▶ Folding function
- ▶ Integrated tilt sensor
- ▶ For measuring yourself and others
- ▶ Measuring range: 50 up to 240 cm



Features

- ▶ Ultrasound measuring technology, temperature-compensated
- ▶ Tilt sensor for correction if measuring instrument is not horizontal, resulting in an error message
- ▶ Low space requirement due to fold function
- ▶ Plastic housing
- ▶ LCD display, 9 mm digit height
- ▶ Acoustic acknowledgment at end of measurement and reading displayed
- ▶ Units in cm or inch
- ▶ Battery operation with 2 x 1.5 V AA batteries (included in scope of delivery)

Products

Ordernumber Variation	Approvable	Division mm	Measuring range cm	Dimensions (mm)
5003.01.001 Height rod, foldable	No	10	50 - 240	345×92×28mm
