



Stand scale with height rod

7831

- ▶ Indication of weight, height and BMI
- ▶ Extra large tread: 320 x 400 mm
- ▶ Integrated height rod
- ▶ Instant-start-function by stepping on the scale
- ▶ Switchable to kg/lb units



Features

- ▶ Integrated height rod
- ▶ Indication of weight, height and BMI
- ▶ Extra large tread: 320 x 400 mm
- ▶ Bright luminous LED display, 15 mm digit height
- ▶ Instant-start-function by stepping on the scale
- ▶ Tare function
- ▶ Mother-child-weighing
- ▶ Easy to clean surface
- ▶ No assembly required - the scale is immediately ready for use
- ▶ Service interface (resp. data cable available)
- ▶ Switchable to kg/lb units
- ▶ Energy-saving auto-off function
- ▶ Mains operation
- ▶ Optional rechargeable battery operation
- ▶ Optional wheels for the easy transport of the scale to the patient

Products

Ordernumber Variation	Approvable	Division g	Maximum load kg	Platform size mm
7831.30.001 Stand scale with height rod	No	100	250	365 x 490