

Foldable IV-Pole

Colour: black

I-H84132

- + simple construction and comfortable handling use for the home care
- + quickly applicable and space-saving storable thanks to its practical travel bag
- +



Product description

Material:	aluminium / plastic
Bottle holder:	2 safety bottle hooks - 0,5 kg / hook
Adjustment:	screw safety height adjustment
Regulating range:	830 - 1.780 mm
Tube diameter:	20 / 15 / 11 mm
Base diameter:	680 mm
Total load capacity:	3 kg
Total weight (net):	appr. 2 kg
Features:	anti-slide stopper base, including carrying bag



Foldable IV-Pole

Colour: black



Application

- + easy and quick to use in emergency and homecare



Technology

- + all plastic parts are light-resisting
- + integrated pull-out guard avoids separation of the adjustable tubes
- + safety bottle holder in accordance with DIN ISO 15375



Economy

- + complete recyclability



Safety

- + enough safe distance in between both tubes once the pole reached its lowest height
- + no risk of injury from the rounded shape of the hook