

IV-Stand Eco-Space

Material: alu / stainless steel

I-S12112

- + effective infusion care by hanging safely IV-bottles and bags
- + mobile use: at the right time at the right place
- + space-saving storage thanks to an intelligent base construction



Product description

Material:	stainless steel
Colour of plastic parts:	blue
Bottle holder:	metal 4 safety bottle hooks - 2 kg / hook
Adjustment:	single handed safety height adjustment
Regulating range:	1.360 - 2.220 mm
Tube diameter:	25 / 18 mm
Length of outer tube:	1.000 mm
Base diameter:	590 mm
Castors:	mono castors Ø 75 mm, according to DIN EN 60601 electrically conductive and lockable
Total / base weight (net):	appr. 3,5 kg / 2,25 kg with drip vessel and holder



IV-Stand Eco-Space

Material: alu / stainless steel



Technology

- + all plastic parts are light-resisting
- + integrated pull-out guard avoids separation of the adjustable tubes
- + safety bottle holder in accordance with DIN ISO 15375

Economy

- + complete recyclability

Safety

- + enough safe distance in between both tubes once the pole reached its lowest height
- + no risk of injury from the rounded shape of the hook

