

MEDICAL DEVICES

JBP Nano Cannula Thin-Walled

A new generation of cannulas have been designed to reduce the damage of hypodermal tissues. Main traits of JBP's Nanocannula are smooth injection and less pain. This thin walled cannula with a wide inner diameter helps insertion with less resistance. The cutting edge design of the cannula hole is micro-polished with a special technology. Each punching needle is designed for less pain and damage to skin tissues.

24 Kits / Box (each kit includes 1 Punching Nanoneedle)



JBP NanoCannula Specification

	JBP NanoCannula Size	Punching Needle Size
30G	Outer Diameter (0.30mm) x 25mm length	27G x 13mm length
27G	Outer Diameter (0.40mm) x 40mm length	26G x 13mm length
27G	Outer Diameter(0.40mm) x 50mm length	26G x 13mm length
25G	Outer Diameter(0.50mm) x 40mm length	23G x 25 mm length
25G	Outer Diameter(0.50mm) x 50mm length	23G x 25mm length
22G	Outer Diameter(0.70mm) x 50mm length	21G x 25mm length
22G	Outer Diameter(0.70mm) x 70mm length	21G x 25mm length

7 different sizes (lengths, diameters and internal calibration) are carefully adapted according to the treated areas, the direction and the depth of injection.

SINGLE USE CE 1293 KFDA KGMP 인증 IFC B000011010485 STERILE EO

 **healthy
medical**
www.hmed.com.hk

Sales Hotline: 8338 3333