

Keyhub 健合®

K5BW Oxygen Concentrator

醫用製氧機 5L/公升

Item code: 003439



- **Compact design with ergonomic handle and 4 wheels**
機身輕巧附有手柄及4輪設計
- **Simple and easy-to-use**
操作簡單，方便易用
- **High-performing and 24 hours continuous operation**
高效能，可連續24小時運作
- **Suitable for home use and medical facilities** 適合家用及診所
- **Comply with CE:Directive 93/42 EEC & ISO13485**
符合歐盟及ISO標準


HEPA Air Intake Filter
配備高效濾網

Efficiency & Energy Saving
高效節能

CE & ISO Standard
符合歐盟及ISO標準



Please contact us at 8338 3333 for details.
請聯絡我們8338 3333 了解產品詳情



Equipped with LED display & knob control
配備LED顯示屏及調節旋轉按鈕

Oxygen Flow Rate	0-5/min
Oxygen Concentration	93%±3% (0-5LPM)
Flow Setting Control	0 to 5 LPM
Maximum Outlet Pressure	8.5psi
Average Sound Level	≤45dB(A)
Power Supply	220 ~ 230V 50HZ
Oxygen Concentration Sensor Alarm	82% or less
Size	305Lx300Wx535H(mm)
Weight	16Kg
Includes	<ol style="list-style-type: none">1. Running time calculation2. Real time oxygen concentration display3. Alarms for power failure, and oxygen flow blockage4. Protection for over-heat and over-load5. Low purity alarm6. Medical nebulization function