





E200 Transport Stretcher

Cogent's transport stretcher provides streamlined design, so it can transport patients efficiently and speedily which has economic value. This transport stretcher is very suitable to be used in Recovery Room and ER. It equipped with four-wheel brake and fifth wheel which offers excellent mobility to provide the best safety to patients. The stretcher also provides the durable quality and 700lbs weight capacity to treat patients with different weights.





Standard Features

- Integrated storage tray with oxygen tank through is made with ABS plastic
- Foot pedal for hydraulic system control on both sides
- Structure with two-section bed board
- Drainage Bag clasp x 2 & Urinary Bag clasp x2
- Foot controlled brake / steer pedal for four wheels in two ends of stretcher
- IV pole inserted hole x5
- Collapsible siderail
- Head-end delivery Push Bar
- Dual pneumatic assisted backrest
- Two stages movable IV pole x 1
- Two stages permanent folding IV pole x 1
- Foot-end assist push bar
- Velcro tape underneath the mattress
- Protective corner bumpers x 4
- Hook / Jack for upright oxygen tank holder
- Retractable steering fifth wheel

Specifications

- Overall Length: 213 cm
- Overall Width: 81/85 cm
- Weight Capacity: 700 lbs (318 kg)
- Height Range (from ground to litter top) :
High: 35" (89 cm) ±1
Low: 22" (54 cm) ±1
- Litter Positioning:
Backrest: 0° - 90°
Trend. / Reverse Trend.: 0-18°
- Patient Surface: 77" x 26" (193 x 64cm) ±2
- Mattress Dimension: 3" (8cm),195 x 66cm±2
- Siderail: 55" x18.7" (141 x48cm)
- Caster Diameter: 8" (20 cm)
- Two years warranty for stretcher under normal use after date of shipping
- Extended warranty available